

AHI1 (H-300): sc-98623

BACKGROUND

AHI1 (abelson helper integration site 1), also known as ORF1, AHI-1, JBTS3 or Jouberin, is a 1,196 amino acid signaling protein that is expressed in the brain, specifically in neurons that give rise to the crossing axons of the corticospinal tract and superior cerebellar peduncles. AHI1 contains seven WD repeats, an SH3 domain, and several SH3-binding sites and is critical for both cerebellar and cortical development. Mutations of AHI1 is associated with Joubert syndrome (JS), an autosomal recessive disorder characterized by hypotonia, ataxia, mental retardation, altered respiratory pattern, abnormal eye movements, and a brain malformation. Considered a novel oncogene, AHI1 is highly deregulated in Chronic myeloid leukemia (CML). Three isoforms exist due to alternative splicing events.

REFERENCES

- Dixon-Salazar, T., Silhavy, J.L., Marsh, S.E., Louie, C.M., Scott, L.C., Gururaj, A., Al-Gazali, L., Al-Tawari, A.A., Kayserili, H., Sztriha, L. and Gleeson, J.G. 2004. Mutations in the AHI1 gene, encoding jouberin, cause Joubert syndrome with cortical polymicrogyria. *Am. J. Hum. Genet.* 75: 979-987.
- Close, J., Game, L., Clark, B., Bergounioux, J., Gerovassili, A. and Thein, S.L. 2004. Genome annotation of a 1.5 Mb region of human chromosome 6q23 encompassing a quantitative trait locus for fetal hemoglobin expression in adults. *BMC Genomics* 5: 33.
- Parisi, M.A., Doherty, D., Eckert, M.L., Shaw, D.W., Ozyurek, H., Aysun, S., Giray, O., Al Swaid, A., Al Shahwan, S., Dohayan, N., Bakhsh, E., Indridason, O.S., Dobyns, W.B., Bennett, C.L., Chance, P.F. and Glass, I.A. 2006. AHI1 mutations cause both retinal dystrophy and renal cystic disease in Joubert syndrome. *J. Med. Genet.* 43: 334-339.
- Parisi, M.A., Doherty, D., Chance, P.F. and Glass, I.A. 2007. Joubert syndrome (and related disorders) (OMIM 213300). *Eur. J. Hum. Genet.* 15: 511-521.
- Kroes, H.Y., van Zon, P.H., Fransen van de Putte, D., Nelen, M.R., Nievelstein, R.J., Wittebol-Post, D., van Nieuwenhuizen, O., Mancini, G.M., van der Knaap, M.S., Kwee, M.L., Maas, S.M., Cobben, J.M., De Nef, J.E., Lindhout, D. and Sinke, R.J. 2008. DNA analysis of AHI1, NPHP1 and CYCLIN D1 in Joubert syndrome patients from the Netherlands. *Eur. J. Med. Genet.* 51: 24-34.

CHROMOSOMAL LOCATION

Genetic locus: AHI1 (human) mapping to 6q23.3; Ahi1 (mouse) mapping to 10 A3.

SOURCE

AHI1 (H-300) is a rabbit polyclonal antibody raised against amino acids 601-897 mapping within an internal region of AHI1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

AHI1 (H-300) is recommended for detection of AHI1 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).


AHI1 (H-300) is also recommended for detection of AHI1 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for AHI1 siRNA (h): sc-72465, AHI1 siRNA (m): sc-72466, AHI1 shRNA Plasmid (h): sc-72465-SH, AHI1 shRNA Plasmid (m): sc-72466-SH, AHI1 shRNA (h) Lentiviral Particles: sc-72465-V and AHI1 shRNA (m) Lentiviral Particles: sc-72466-V.

Molecular Weight of AHI1: 136 kDa.

Positive Controls: Jurkat whole cell lysate: sc-2204, HEK293 whole cell lysate: sc-45136 or NTERA-2 cl.D1 whole cell lysate: sc-364181.

DATA



AHI1 (H-300): sc-98623. Western blot analysis of AHI1 expression in Jurkat (A), HEK293 (B), NTERA-2 cl.D1 (C), IMR-32 (D) and SH-SY5Y (E) whole cell lysates and mouse testis tissue extract (F).

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **AHI1 (B-2): sc-515382**, our highly recommended monoclonal alternative to AHI1 (H-300).