



## Scabrous (D-300): sc-98848

### BACKGROUND

*Drosophila melanogaster* is a proven and effective model for studying developmental and cellular processes common to higher eukaryotes. Approximately 13,600 genes have been elucidated from more than 120 megabases of euchromatin, and they are organized among the chromosomes 2, 3, 4, X and Y, with the Y chromosome being predominately heterochromatic. *Drosophila* genes can be categorized based on the type of protein for which they encode and are represented by six major classifications, which include intracellular signaling proteins, transmembrane proteins, RNA binding proteins, secreted factors, transcription regulators (basic helix-loop-helix, homeodomain-containing, zinc finger-containing, and chromatin-associated) or other functional proteins. *Drosophila* scabrous (*sca*) gene encodes a secreted dimeric glycoprotein with putative coiled-coil domains N-terminally and a C-terminal region related to the blood clot protein Fibrinogen. Scabrous influences ommatidial rotation during eye development by suppressing nemo activity. Homozygous *sca* mutants have extra bristle organs and rough eyes.

### REFERENCES

1. Mlodzik, M., Baker, N.E. and Rubin, G.M. 1990. Isolation and expression of Scabrous, a gene regulating neurogenesis in *Drosophila*. *Genes Dev.* 4: 1848-1861.
2. Ellis, M.C., Weber, U., Wiersdorff, V. and Mlodzik, M. 1994. Confrontation of Scabrous expressing and non-expressing cells is essential for normal ommatidial spacing in the *Drosophila* eye. *Development* 120: 1959-1969.
3. Lai, C., Lyman, R.F., Long, A.D., Langley, C.H. and Mackay, T.F. 1994. Naturally occurring variation in bristle number and DNA polymorphisms at the Scabrous locus of *Drosophila melanogaster*. *Science* 266: 1697-1702.
4. Hu, X., Lee, E.C. and Baker, N.E. 1995. Molecular analysis of Scabrous mutant alleles from *Drosophila melanogaster* indicates a secreted protein with two functional domains. *Genetics* 141: 607-617.
5. Lee, E.C., Hu, X., Yu, S.Y. and Baker, N.E. 1996. The Scabrous gene encodes a secreted glycoprotein dimer and regulates proneural development in *Drosophila* eyes. *Mol. Cell Biol.* 16: 1179-1188.
6. Lee, E.C., Yu, S.Y., Hu, X., Mlodzik, M. and Baker, N.E. 1998. Functional analysis of the Fibrinogen-related scabrous gene from *Drosophila melanogaster* identifies potential effector and stimulatory protein domains. *Genetics* 150: 663-673.
7. Adams, M.D., Celniker, S.E., Holt, R.A., Evans, C.A., Gocayne, J.D., Amanatides, P.G., Scherer, S.E., Li, P.W., Hoskins, R.A., Galle, R.F., George, R.A., Lewis, S.E., Richards, S., Ashburner, M., Henderson, S.N., et al. 2000. The genome sequence of *Drosophila melanogaster*. *Science* 287: 2185-2195.
8. Renaud, O. and Simpson, P. 2001. Scabrous modifies epithelial cell adhesion and extends the range of lateral signalling during development of the spaced bristle pattern in *Drosophila*. *Dev. Biol.* 240: 361-376.
9. Chou, Y.H. and Chien, C.T. 2002. Scabrous controls ommatidial rotation in the *Drosophila* compound eye. *Dev. Cell* 3: 839-850.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### SOURCE

Scabrous (D-300) is a rabbit polyclonal antibody raised against amino acids 468-767 mapping near the C-terminus of Scabrous of *Drosophila melanogaster* origin.

### PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

### APPLICATIONS

Scabrous (D-300) is recommended for detection of Scabrous of *Drosophila melanogaster* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

### RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

### STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.