

RAI1 (H-270): sc-98883

BACKGROUND

Retinoic acid induced 1 (RAI1) is a 1,906 amino acid protein containing an N-terminal polyglutamine stretch that is expressed in most tissues, with highest expression in neuronal tissues. RAI1 functions as a transcriptional regulator and is important for embryonic and postnatal developments. Heterozygous deletions of the RAI1 gene are associated with Smith-Magenis syndrome (SMS), a mental retardation syndrome with behavioral, neurological and skeletal anomalies. Individuals affected with SMS usually display self-injurious behaviors, sleep disturbance, developmental delay and reduced motor and cognitive skills. RAI1 haploinsufficiency is specifically responsible for the obesity and craniofacial symptoms of SMS. RAI1 mutations have also been implicated in schizophrenia and spinocerebellar ataxia type 2.

REFERENCES

1. Rebohle, E., Reimer, W., Berg, E., Buchecker, H., Gewalt, I., Klingbeil, M., Köhler, E., Lindemann, R., Meinert, G., Mencik, M., Nack, P., Reum, P.J. and Schneider, D. 1977. Results of additional studies on the performance-diagnostic screening test. *Z. Gesamte Hyg.* 23: 896-899.
2. Hayes, S., Turecki, G., Brisebois, K., Lopes-Cendes, I., Gaspar, C., Riess, O., Ranum, L.P., Pulst, S.M. and Rouleau, G.A. 2000. CAG repeat length in RAI1 is associated with age at onset variability in spinocerebellar ataxia type 2 (SCA2). *Hum. Mol. Genet.* 9: 1753-1758.
3. Slager, R.E., Newton, T.L., Vlangos, C.N., Finucane, B. and Elsea, S.H. 2003. Mutations in RAI1 associated with Smith-Magenis syndrome. *Nat. Genet.* 33: 466-468.
4. Toulouse, A., Rochefort, D., Rousset, J., Joobert, R. and Rouleau, G.A. 2003. Molecular cloning and characterization of human RAI1, a gene associated with schizophrenia. *Genomics* 82: 162-171.
5. Bi, W., Saifi, G.M., Shaw, C.J., Walz, K., Fonseca, P., Wilson, M., Potocki, L. and Lupski, J.R. 2004. Mutations of RAI1, a PHD-containing protein, in nondeletion patients with Smith-Magenis syndrome. *Hum. Genet.* 115: 515-524.
6. Bi, W., Ohyama, T., Nakamura, H., Yan, J., Visvanathan, J., Justice, M.J. and Lupski, J.R. 2005. Inactivation of RAI1 in mice recapitulates phenotypes observed in chromosome engineered mouse models for Smith-Magenis syndrome. *Hum. Mol. Genet.* 14: 983-995.
7. Girirajan, S., Elsas, L.J., Devriendt, K. and Elsea, S.H. 2005. RAI1 variations in Smith-Magenis syndrome patients without 17p11.2 deletions. *J. Med. Genet.* 42: 820-828.
8. Vlangos, C.N., Wilson, M., Blancato, J., Smith, A.C. and Elsea, S.H. 2005. Diagnostic FISH probes for del(17)(p11.2;p11.2) associated with Smith-Magenis syndrome should contain the RAI1 gene. *Am. J. Med. Genet. A* 132: 278-282.

CHROMOSOMAL LOCATION

Genetic locus: RAI1 (human) mapping to 17p11.2; Rai1 (mouse) mapping to 11 B1.3.

SOURCE

RAI1 (H-270) is a rabbit polyclonal antibody raised against amino acids 6-275 mapping near the N-terminus of RAI1 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

RAI1 (H-270) is recommended for detection of RAI1 isoforms 1, 2, 3 and 4 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for RAI1 siRNA (h): sc-61438, RAI1 siRNA (m): sc-61439, RAI1 shRNA Plasmid (h): sc-61438-SH, RAI1 shRNA Plasmid (m): sc-61439-SH, RAI1 shRNA (h) Lentiviral Particles: sc-61438-V and RAI1 shRNA (m) Lentiviral Particles: sc-61439-V.

Molecular Weight of RAI1: 203 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **RAI1 (D-11): sc-365065**, our highly recommended monoclonal alternative to RAI1 (H-270).