

# DOT1 (Y-290): sc-98944



The Power to Question

## BACKGROUND

Yeast disruptor of telomeric silencing-1 (DOT1) functions in gene silencing at the pachytene checkpoint during the meiotic cell cycle. This owes to the fact that DOT1 methylates Histone H3 at Lysine 79. Methylation of histone proteins affects chromatin structure to regulate gene expression, in this case, resulting in silencing of genes located near chromosome telomeres. This methylation limits silencing to discrete loci by preventing the binding of Sir proteins along the genome, and as these interacting proteins are conserved across species, similar mechanisms likely exist in other eukaryotes.

## REFERENCES

- Gottschling, D.E. 1998. Identification of high-copy disruptors of telomeric silencing in *Saccharomyces cerevisiae*. *Genetics* 150: 613-632.
- San-Segundo, P.A. and Roeder, G.S. 2000. Role for the silencing protein DOT1 in meiotic checkpoint control. *Mol. Biol. Cell* 11: 3601-3615.
- Lacoste, N., Utley, R.T., Hunter, J.M., Poirier, G.G. and Cote, J. 2002. Disruptor of telomeric silencing-1 is a chromatin-specific Histone H3 methyltransferase. *J. Biol. Chem.* 277: 30421-30424.
- van Leeuwen, F., Gafken, P.R. and Gottschling, D.E. 2002. DOT1p modulates silencing in yeast by methylation of the nucleosome core. *Cell* 109: 745-756.
- Ng, H.H., Feng, Q., Wang, H., Erdjument-Bromage, H., Tempst, P., Zhang, Y. and Struhl, K. 2002. Lysine methylation within the globular domain of Histone H3 by DOT1 is important for telomeric silencing and Sir protein association. *Genes Dev.* 16: 1518-1527.
- Ng, H.H., Xu, R.M., Zhang, Y. and Struhl, K. 2002. Ubiquitination of Histone H2B by Rad6 is required for efficient DOT1-mediated methylation of Histone H3 Lysine 79. *J. Biol. Chem.* 277: 34655-34657.
- Krogan, N.J., Dover, J., Wood, A., Schneider, J., Heidt, J., Boateng, M.A., Dean, K., Ryan, O.W., Golshani, A., Johnston, M., Greenblatt, J.F. and Shilatifard, A. 2003. The Paf1 complex is required for Histone H3 methylation by COMPASS and DOT1p: linking transcriptional elongation to histone methylation. *Mol. Cell* 11: 721-729.
- Zhang, W., Hayashizaki, Y. and Kone, B.C. 2004. Structure and regulation of the mDot1 gene, a mouse histone H3 methyltransferase. *Biochem J.* 377: 641-651.

## SOURCE

DOT1 (Y-290) is a rabbit polyclonal antibody raised against amino acids 293-582 mapping at the C-terminus of DOT1 of *Saccharomyces cerevisiae* origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## APPLICATIONS

DOT1 (Y-290) is recommended for detection of DOT1 of *Saccharomyces cerevisiae* origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.