

# NR3A (H-190): sc-98986

## BACKGROUND

NR3A is a subunit of the N-methyl-D-aspartate (NMDA) receptors, which belong to the superfamily of glutamate-regulated ion channels and function in pathological and physiological processes in the central nervous system. NR3A is a multi-pass membrane protein that is expressed in fetal brain and is mediated by glycine. It may be involved in the development of dendritic spines and in the PPP2CB-NMDAR mediated signaling mechanism. NR3A forms a heteromeric channel composed of a  $\zeta$  subunit (GRIN1), an  $\epsilon$  subunit (GRIN2A, GRIN2B, GRIN2C or GRIN2D) and a third subunit (GRIN3A or GRIN3B). The NR3A protein is enriched in post-synaptic plasma membrane and post-synaptic densities and requires the presence of GRIN1 to be targeted at the plasma membrane. The NR3A subunit displays greater than 90% sequence homology to the corresponding subunit in rat.

## REFERENCES

1. Das, S., Sasaki, Y.F., Rothe, T., Premkumar, L.S., Takasu, M., Crandall, J.E., Dikkes, P., Conner, D.A., Rayudu, P.V., Cheung, W., Chen, H.S., Lipton, S.A. and Nakanishi, N. 1998. Increased NMDA current and spine density in mice lacking the NMDA receptor subunit NR3A. *Nature* 393: 377-381.
2. Chan, S.F. and Sucher, N.J. 2001. An NMDA receptor signaling complex with protein phosphatase 2A. *J. Neurosci.* 21: 7985-7992.
3. Andersson, O., Stenqvist, A., Attersand, A. and von Euler, G. 2001. Nucleotide sequence, genomic organization, and chromosomal localization of genes encoding the human NMDA receptor subunits NR3A and NR3B. *Genomics* 78: 178-184.
4. Chatterton, J.E., Awobuluyi, M., Premkumar, L.S., Takahashi, H., Talantova, M., Shin, Y., Cui, J., Tu, S., Sevarino, K.A., Nakanishi, N., Tong, G., Lipton, S.A. and Zhang, D. 2002. Excitatory glycine receptors containing the NR3 family of NMDA receptor subunits. *Nature* 415: 793-798.
5. Eriksson, M., Nilsson, A., Froelich-Fabre, S., Akesson, E., Dunker, J., Seiger, A., Folkesson, R., Benedikz, E. and Sundström, E. 2002. Cloning and expression of the human N-methyl-D-aspartate receptor subunit NR3A. *Neurosci. Lett.* 321: 177-181.
6. Nagase, T., Kikuno, R. and Ohara, O. 2002. Prediction of the coding sequences of unidentified human genes. XXII. The complete sequences of 50 new cDNA clones which code for large proteins. *DNA Res.* 8: 319-327.
7. Mueller, H.T. and Meador-Woodruff, J.H. 2003. Expression of the NR3A subunit of the NMDA receptor in human fetal brain. *Ann. N.Y. Acad. Sci.* 1003: 448-451.

## CHROMOSOMAL LOCATION

Genetic locus: GRIN3A (human) mapping to 9q31.1; Grin3a (mouse) mapping to 4 B1.

## SOURCE

NR3A (H-190) is a rabbit polyclonal antibody raised against amino acids 201-390 mapping within an N-terminal extracellular domain of NR3A of human origin.

## PRODUCT

Each vial contains 200  $\mu$ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

NR3A (H-190) is recommended for detection of NR3A of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2  $\mu$ g per 100-500  $\mu$ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

NR3A (H-190) is also recommended for detection of NR3A in additional species, including equine, canine and bovine.

Suitable for use as control antibody for NR3A siRNA (h): sc-61229, NR3A siRNA (m): sc-61230, NR3A shRNA Plasmid (h): sc-61229-SH, NR3A shRNA Plasmid (m): sc-61230-SH, NR3A shRNA (h) Lentiviral Particles: sc-61229-V and NR3A shRNA (m) Lentiviral Particles: sc-61230-V.

Molecular Weight of NR3A: 135 kDa.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) or our catalog for detailed protocols and support products.