

BBS7 (H-300): sc-99094

BACKGROUND

Bardet-Biedl syndrome (BBS) is a pleiotropic genetic disorder characterized by obesity, photoreceptor degeneration, polydactyly, hypogonadism, renal abnormalities and developmental delay. BBS patients also have an increased risk of developing diabetes, hypertension and congenital heart defects. BBS is a heterogeneous disorder; BBS genes map to eight genetic loci and encode eight proteins, BBS1-BBS8. Five BBS genes encode basal body or cilia proteins, suggesting that BBS is a ciliary dysfunction disorder. The BBS2 gene contains two overlapping genes: BBS2L1 and BBS2L2. BBSL1 was renamed BBS7, whereas BBS2L2 independently functions as BBS1. BBS7 contains 672 amino acids and is expressed at low to moderate levels in most human tissues.

CHROMOSOMAL LOCATION

Genetic locus: BBS7 (human) mapping to 4q27; Bbs7 (mouse) mapping to 3 B.

SOURCE

BBS7 (H-300) is a rabbit polyclonal antibody raised against amino acids 416-715 mapping at the C-terminus of BBS7 of human origin.

PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin. Also available as TransCruz reagent for Gel Supershift and ChIP applications, sc-99094 X, 200 µg/0.1 ml.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

BBS7 (H-300) is recommended for detection of BBS7 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

BBS7 (H-300) is also recommended for detection of BBS7 in additional species, including equine, canine, bovine and porcine.

Suitable for use as control antibody for BBS7 siRNA (h): sc-60259, BBS7 siRNA (m): sc-60260, BBS7 shRNA Plasmid (h): sc-60259-SH, BBS7 shRNA Plasmid (m): sc-60260-SH, BBS7 shRNA (h) Lentiviral Particles: sc-60259-V and BBS7 shRNA (m) Lentiviral Particles: sc-60260-V.

BBS7 (H-300) X TransCruz antibody is recommended for Gel Supershift and ChIP applications.


Molecular Weight of BBS7: 80 kDa.

Positive Controls: BBS7 (m): 293T Lysate: sc-118689.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

DATA



BBS7 (H-300): sc-99094. Western blot analysis of BBS7 expression in non-transfected: sc-117752 (A) and mouse BBS7 transfected: sc-118689 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.



Try **BBS7 (E-8): sc-390403**, our highly recommended monoclonal alternative to BBS7 (H-300).