

# CLCA2/5 (H-60): sc-99223

## BACKGROUND

The calcium-activated chloride channel (CLCA) protein family, which includes the human homologs CLCA1 and CLCA2, display distinct tissue distribution patterns. CLCA1 is expressed as a precursor protein that is processed into two cell surface associated subunits and a group of proteins. CLCA1 is upregulated by interleukin-9 and regulates the expression of mucins. CLCA1 may provide a therapeutic target to control mucus overproduction in airway disease patients with cystic fibrosis. CLCA2 expression is downregulated in breast cancer and, therefore, is thought to act as a tumor suppressor in normal cells. CLCA3 is a structurally divergent member of the CLCA family that does not function as a channel protein. CLCA4 is a CLCA member that is expressed in human rectal mucosa. CLCA5 shows strong expression in eye and spleen, and CLCA6 is primarily expressed in intestine and stomach.

## REFERENCES

- Gandhi, R., Elble, R.C., Gruber, A.D., Schreur, K.D., Ji, H.L., Fuller, C.M. and Pauli, B.U. 1998. Molecular and functional characterization of a calcium-sensitive chloride channel from mouse lung. *J. Biol. Chem.* 273: 32096-32101.
- Gruber, A.D., Elble, R.C., Ji, H.L., Schreur, K.D., Fuller, C.M. and Pauli, B.U. 1998. Genomic cloning, molecular characterization, and functional analysis of human CLCA1, the first human member of the family of Ca<sup>2+</sup>-activated Cl<sup>-</sup> channel proteins. *Genomics* 54: 200-214.
- Gruber, A.D., Schreur, K.D., Ji, H.L., Fuller, C.M. and Pauli, B.U. 1999. Molecular cloning and transmembrane structure of hCLCA2 from human lung, trachea, and mammary gland. *Am. J. Physiol.* 276: 1261-1270.
- Hauber, H.P., Manoukian, J.J., Nguyen, L.H., Sobol, S.E., Levitt, R.C., Holroyd, K.J., McElvaney, N.G., Griffin, S. and Hamid, Q. 2003. Increased expression of interleukin-9, interleukin-9 receptor, and the calcium-activated chloride channel hCLCA1 in the upper airways of patients with cystic fibrosis. *Laryngoscope* 113: 1037-1042.
- Beckley, J.R., Pauli, B.U. and Elble, R.C. 2004. Re-expression of detachment-inducible chloride channel mCLCA5 suppresses growth of metastatic breast cancer cells. *J. Biol. Chem.* 279: 41634-41641.
- Hauber, H.P., Tscopoulos, A., Wallaert, B., Griffin, S., McElvaney, N.G., Daigneault, P., Mueller, Z., Olivenstein, R., Holroyd, K.J., Levitt, R.C. and Hamid, Q. 2004. Expression of hCLCA1 in cystic fibrosis lungs is associated with mucus overproduction. *Eur. Respir. J.* 23: 846-850.
- Li, X., Cowell, J.K. and Sossey-Alaoui, K. 2004. CLCA2 tumour suppressor gene in 1p31 is epigenetically regulated in breast cancer. *Oncogene* 23: 1474-1480.

## CHROMOSOMAL LOCATION

Genetic locus: CLCA2 (human) mapping to 1p22.3; Clca5 (mouse) mapping to 3 H2.

## SOURCE

CLCA2/5 (H-60) is a rabbit polyclonal antibody raised against amino acids 601-660 mapping within an internal region of CLCA2 of human origin.

## PRODUCT

Each vial contains 200 µg IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

CLCA2/5 (H-60) is recommended for detection of CLCA2 of human origin and, to a lesser extent, Clca5 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CLCA2/5 (H-60) is also recommended for detection of CLCA2 in additional species, including porcine.

Suitable for use as control antibody for CLCA2 siRNA (h): sc-60394, Clca5 siRNA (m): sc-60395, CLCA2 shRNA Plasmid (h): sc-60394-SH, Clca5 shRNA Plasmid (m): sc-60395-SH, CLCA2 shRNA (h) Lentiviral Particles: sc-60394-V and Clca5 shRNA (m) Lentiviral Particles: sc-60395-V.

Molecular Weight of CLCA2/5 precursor: 120 kDa.

Molecular Weight of CLCA2/5 subunits: 86/34 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200, HUV-EC-C + VEGF cell lysate: sc-24709 or MCF7 whole cell lysate: sc-2206.

## RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use goat anti-rabbit IgG-HRP: sc-2004 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible goat anti-rabbit IgG-HRP: sc-2030 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use goat anti-rabbit IgG-FITC: sc-2012 (dilution range: 1:100-1:400) or goat anti-rabbit IgG-TR: sc-2780 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.



Try **CLCA2 (1D5): sc-517015** or **Clca5 (A-5): sc-515682**, our highly recommended monoclonal alternatives to CLCA2/5 (H-60).