

IDI1 (XY-7): sc-100550

BACKGROUND

IDI1 (isopentenyl-diphosphate δ isomerase 1), also known as IPP1 or IPP11, is a 227 amino acid member of the IPP isomerase type I family and is involved in cholesterol biosynthesis. Localized to the peroxisome, IDI1 catalytically converts isopentenyl diphosphate (IPP) to its electrophilic isomer, dimethylallyl diphosphate (DMAPP). Specifically, IDI1 uses magnesium as a cofactor to catalyze the 1,3-allylic rearrangement of IPP, thus creating DMAPP, a substrate for subsequent reactions that synthesize farnesyl diphosphate and, ultimately, cholesterol. Defects in the gene encoding IDI1 may be associated with peroxisomal deficiency diseases, such as Zellweger syndrome, a congenital disorder caused by a reduction in the number of peroxisomes. Individuals affected with this disorder generally exhibit lack of muscle tone, an enlarged liver, mental retardation and, in some cases, death.

REFERENCES

- Xuan, J.W., Kowalski, J., Chambers, A.F. and Denhardt, D.T. 1994. A human promyelocyte mRNA transiently induced by TPA is homologous to yeast IPP isomerase. *Genomics* 20: 129-131.
- Krisans, S.K., Ericsson, J., Edwards, P.A. and Keller, G.A. 1994. Farnesyl-diphosphate synthase is localized in peroxisomes. *J. Biol. Chem.* 269: 14165-14169.
- Hahn, F.M., Xuan, J.W., Chambers, A.F. and Poulter, C.D. 1996. Human isopentenyl-diphosphate: dimethylallyl diphosphate isomerase: overproduction, purification, and characterization. *Arch. Biochem. Biophys.* 332: 30-34.
- Paton, V.G., Shackelford, J.E. and Krisans, S.K. 1997. Cloning and subcellular localization of hamster and rat isopentenyl-diphosphate dimethylallyl diphosphate isomerase. A PTS1 motif targets the enzyme to peroxisomes. *J. Biol. Chem.* 272: 18945-18950.
- Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604055. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Breitling, R., Laubner, D., Clizbe, D., Adamski, J. and Krisans, S.K. 2003. Isopentenyl-diphosphate isomerases in human and mouse: evolutionary analysis of a mammalian gene duplication. *J. Mol. Evol.* 57: 282-291.
- Clizbe, D.B., Owens, M.L., Masuda, K.R., Shackelford, J.E. and Krisans, S.K. 2007. IDI2, a second isopentenyl-diphosphate isomerase in mammals. *J. Biol. Chem.* 282: 6668-6676.
- Zheng, W., Sun, F., Bartlam, M., Li, X., Li, R. and Rao, Z. 2007. The crystal structure of human isopentenyl-diphosphate isomerase at 1.7 Å resolution reveals its catalytic mechanism in isoprenoid biosynthesis. *J. Mol. Biol.* 366: 1447-1458.

CHROMOSOMAL LOCATION

Genetic locus: IDI1 (human) mapping to 10p15.3.

SOURCE

IDI1 (XY-7) is a mouse monoclonal antibody raised against recombinant IDI1 of human origin.

PRODUCT

Each vial contains 100 μ g IgG₁ kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

IDI1 (XY-7) is recommended for detection of IDI1 of human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)] and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for IDI1 siRNA (h): sc-90838, IDI1 shRNA Plasmid (h): sc-90838-SH and IDI1 shRNA (h) Lentiviral Particles: sc-90838-V.


Molecular Weight of IDI1: 26 kDa.

Positive Controls: HeLa whole cell lysate: sc-2200.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml).

DATA



IDI1 (XY-7): sc-100550. Western blot analysis of IDI1 expression in HeLa whole cell lysate.

SELECT PRODUCT CITATIONS

- Brisdelli, F., Di Francesco, L., Giorgi, A., Lizzi, A.R., Luzi, C., Mignogna, G., Bozzi, A. and Schininà, M.E. 2019. Proteomic analysis of quercetin-treated K562 cells. *Int. J. Mol. Sci.* 21 pii: E32.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.