

FKBP9 siRNA (m): sc-105356

BACKGROUND

The immunophilins are a highly conserved family of *cis-trans* peptidyl-prolyl isomerases that bind to and mediate the effects of immunosuppressive drugs, such as cyclosporin, FK506 and rapamycin. Immunophilins have also been implicated in protein folding and trafficking within the endoplasmic reticulum. FKBP9 (FK506-binding protein 9), also known as Rotamase, 63 kDa FK506-binding protein and peptidyl-prolyl *cis-trans* isomerase FKBP9, is a 570 amino acid endoplasmic reticular protein that belongs to the FKBP-type PPIase family, a group of proteins known to catalyze the folding of proline-containing polypeptides. FKBP9 is expressed at high levels in mouse lung, kidney, heart and muscle. The gene encoding FKBP9 maps to human chromosome 7, which houses over 1,000 genes and comprises nearly 5% of the human genome. Defects in some of the genes localized to chromosome 7 have been linked to osteogenesis imperfecta, Williams-Beuren syndrome, Pendred syndrome, lissencephaly, citrullinemia and Shwachman-Diamond syndrome.

REFERENCES

1. Gold, B.G. 1997. FK506 and the role of immunophilins in nerve regeneration. *Mol. Neurobiol.* 15: 285-306.
2. Shadidy, M., Caubit, X., Olsen, R., Seternes, O.M., Moens, U. and Krauss, S. 1999. Biochemical analysis of mouse FKBP60, a novel member of the FKPB family. *Biochim. Biophys. Acta* 1446: 295-307.
3. Jo, D., Lyu, M.S., Cho, E.G., Park, D., Kozak, C.A. and Kim, M.G. 2001. Identification and genetic mapping of the mouse Fkbp9 gene encoding a new member of FK506-binding protein family. *Mol. Cells* 12: 272-275.
4. Gilbert, F. 2002. Chromosome 7. *Genet. Test.* 6: 141-161.
5. Hillier, L.W., Fulton, R.S., Fulton, L.A., Graves, T.A., Pepin, K.H., Wagner-McPherson, C., Layman, D., Maas, J., Jaeger, S., Walker, R., Wylie, K., Sekhon, M., Becker, M.C., O'Laughlin, M.D., et al. 2003. The DNA sequence of human chromosome 7. *Nature* 424: 157-164.
6. Bhuiyan, Z.A., Hamdan, M.A., Shamsi, E.T., Postma, A.V., Mannens, M.M., Wilde, A.A. and Al-Gazali, L. 2007. A novel early onset lethal form of catecholaminergic polymorphic ventricular tachycardia maps to chromosome 7p14-p22. *J. Cardiovasc. Electrophysiol.* 18: 1060-1066.
7. Kang, C.B., Hong, Y., Dhe-Paganon, S. and Yoon, H.S. 2008. FKBP family proteins: immunophilins with versatile biological functions. *Neurosignals* 16: 318-325.

CHROMOSOMAL LOCATION

Genetic locus: Fkbp9 (mouse) mapping to 6 B3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

FKBP9 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see FKBP9 shRNA Plasmid (m): sc-105356-SH and FKBP9 shRNA (m) Lentiviral Particles: sc-105356-V as alternate gene silencing products.

For independent verification of FKBP9 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-105356A, sc-105356B and sc-105356C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

FKBP9 siRNA (m) is recommended for the inhibition of FKBP9 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor FKBP9 gene expression knockdown using RT-PCR Primer: FKBP9 (m)-PR: sc-105356-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.