

Histone cluster 1 H2BJ siRNA (h): sc-105480

BACKGROUND

Eukaryotic histones are basic and water soluble nuclear proteins that form hetero-octameric nucleosome particles by wrapping 146 base pairs of DNA in a left-handed super-helical turn, sequentially, to form chromosomal fibers. Two molecules of each of the four core histones (Histone H2A, H2B, H3, and H4) form the octamer, which consists of two H2A-H2B dimers and two H3-H4 dimers that are nearly symmetrical by tertiary structure. The Histone H2B core protein is accompanied by several histone cluster proteins, including Histone cluster 1 H2BK, Histone cluster 1 H2BI and Histone cluster 1 H2BJ, all of which help to facilitate chromatin assembly and may exhibit antibacterial activity. Over 80% of nucleosomes contain the linker Histone H1, derived from an intronless gene, that interacts with linker DNA between nucleosomes and mediates compaction into higher order chromatin. Histones are subject to posttranslational modifications, including methylation, citrullination, acetylation, phosphorylation, sumoylation, ubiquitination and ADP-ribosylation.

REFERENCES

- Zhong, R., et al. 1983. The primary structure and expression of four cloned human histone genes. *Nucleic Acids Res.* 11: 7409-7425.
- Albig, W., et al. 1999. The human H2A and H2B Histone gene complement. *Biol. Chem.* 380: 7-18.
- Kim, H.S., et al. 2002. Endotoxin-neutralizing antimicrobial proteins of the human placenta. *J. Immunol.* 168: 2356-2364.
- Cheung, W.L., et al. 2003. Apoptotic phosphorylation of Histone H2B is mediated by mammalian sterile twenty kinase. *Cell* 113: 507-517.
- Tollin, M., et al. 2003. Antimicrobial peptides in the first line defence of human colon mucosa. *Peptides* 24: 523-530.
- Zhu, B., et al. 2005. Monoubiquitination of human Histone H2B: the factors involved and their roles in HOX gene regulation. *Mol. Cell* 20: 601-611.
- Pavri, R., et al. 2006. Histone H2B monoubiquitination functions cooperatively with FACT to regulate elongation by RNA polymerase II. *Cell* 125: 703-717.
- Bonenfant, D., et al. 2006. Characterization of Histone H2A and H2B variants and their post-translational modifications by mass spectrometry. *Mol. Cell. Proteomics* 5: 541-552.
- Beck, H.C., et al. 2006. Quantitative proteomic analysis of post-translational modifications of human histones. *Mol. Cell. Proteomics* 5: 1314-1325.

CHROMOSOMAL LOCATION

Genetic locus: HIST1H2BJ (human) mapping to 6p22.1.

PRODUCT

Histone cluster 1 H2BJ siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Histone cluster 1 H2BJ shRNA Plasmid (h): sc-105480-SH and Histone cluster 1 H2BJ shRNA (h) Lentiviral Particles: sc-105480-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCL, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Histone cluster 1 H2BJ siRNA (h) is recommended for the inhibition of Histone cluster 1 H2BJ expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.