

# Histone cluster 1 H3E siRNA (h): sc-105490

## BACKGROUND

Eukaryotic histones are basic and water soluble nuclear proteins that form hetero-octameric nucleosome particles by wrapping 146 base pairs of DNA in a left-handed super-helical turn sequentially to form chromosomal fiber. Two molecules of each of the four core Histones (H2A, H2B, H3, and H4) form the octamer, which is comprised of two H2A-H2B dimers and two H3-H4 dimers, forming two nearly symmetrical halves by tertiary structure. Histones are subject to posttranslational modification by enzymes primarily on their N-terminal tails, but also in their globular domains. Histone cluster 1 H3E is a member of the H3 Histone family and is encoded by a gene located in the large histone gene cluster on chromosome 6p22.2. Histone cluster 1 H3E undergoes posttranslational modifications, including the phosphorylation of Ser 10.

## REFERENCES

1. Albig, W., Kardalinos, E., Drabent, B., Zimmer, A. and Doenecke, D. 1991. Isolation and characterization of two human H1 Histone genes within clusters of core histone genes. *Genomics* 10: 940-948.
2. Albig, W., Kioschis, P., Poustka, A., Meergans, K. and Doenecke, D. 1997. Human histone gene organization: nonregular arrangement within a large cluster. *Genomics* 40: 314-322.
3. Albig, W. and Doenecke, D. 1997. The human histone gene cluster at the D6S105 locus. *Hum. Genet.* 101: 284-294.
4. Goto, H., Tomono, Y., Ajiro, K., Kosako, H., Fujita, M., Sakurai, M., Okawa, K., Iwamatsu, A., Okigaki, T., Takahashi, T. and Inagaki, M. 1999. Identification of a novel phosphorylation site on Histone H3 coupled with mitotic chromosome condensation. *J. Biol. Chem.* 274: 25543-25549.
5. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602813. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Wang, H., Zhai, L., Xu, J., Joo, H.Y., Jackson, S., Erdjument-Bromage, H., Tempst, P., Xiong, Y. and Zhang, Y. 2006. Histone H3 and H4 ubiquitylation by the CUL4-DDB-ROC1 ubiquitin ligase facilitates cellular response to DNA damage. *Mol. Cell* 22: 383-394.
7. Beck, H.C., Nielsen, E.C., Matthiesen, R., Jensen, L.H., Sehested, M., Finn, P., Grauslund, M., Hansen, A.M. and Jensen, O.N. 2006. Quantitative proteomic analysis of post-translational modifications of human histones. *Mol. Cell. Proteomics* 5: 1314-1325.
8. Morris, S.A., Rao, B., Garcia, B.A., Hake, S.B., Diaz, R.L., Shabanowitz, J., Hunt, D.F., Allis, C.D., Lieb, J.D. and Strahl, B.D. 2007. Identification of Histone H3 lysine 36 acetylation as a highly conserved histone modification. *J. Biol. Chem.* 282: 7632-7640.
9. Garcia, B.A., Hake, S.B., Diaz, R.L., Kauer, M., Morris, S.A., Recht, J., Shabanowitz, J., Mishra, N., Strahl, B.D., Allis, C.D. and Hunt, D.F. 2007. Organismal differences in post-translational modifications in Histones H3 and H4. *J. Biol. Chem.* 282: 7641-7655.

## CHROMOSOMAL LOCATION

Genetic locus: HIST1H3E (human) mapping to 6p22.2.

## PRODUCT

Histone cluster 1 H3E siRNA (h) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Histone cluster 1 H3E shRNA Plasmid (h): sc-105490-SH and Histone cluster 1 H3E shRNA (h) Lentiviral Particles: sc-105490-V as alternate gene silencing products.

For independent verification of Histone cluster 1 H3E (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-105490A and sc-105490B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Histone cluster 1 H3E siRNA (h) is recommended for the inhibition of Histone cluster 1 H3E expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.