



MCT6 siRNA (h): sc-106212

BACKGROUND

Monocarboxylates, such as lactate and pyruvate, play an integral role in cellular metabolism. Lactic acid is produced in large quantities as a result of glycolysis, which provides the majority of ATP to cells under normal physiological conditions. However, accumulation of lactic acid leads to a decrease in intracellular pH and cessation of glycolysis. In order for glycolysis to continue at a high rate, lactic acid must be transported out of the cell. This transport process is carried out by a family of monocarboxylate transporters (MCTs), which function as proton symports and are stereoselective for L-lactate. The MCT family consists of at least 8 members, MCT1-8, which contain between 10-12 transmembrane-helical (TM) domains, with the amino- and carboxy-termini located in the cytoplasm. MCT6 is highly expressed in the kidneys and is thought to have a specificity for bumetanide, a loop diuretic involved in the treatment of edema. The high substrate specificity of MCT6 suggests a possible role in therapeutic drug transport and trafficking across the plasma membrane.

REFERENCES

- Price, N.T., Jackson, V.N. and Halestrap, A.P. 1998. Cloning and sequencing of four new mammalian monocarboxylate transporter (MCT) homologues confirms the existence of a transporter family with an ancient past. *Biochem. J.* 329: 321-328.
- Hadjigapriou, C., Schmidt, L., Dudeja, P.K., Layden, T.J. and Ramaswamy, K. 2000. Mechanism(s) of butyrate transport in Caco-2 cells: role of monocarboxylate transporter 1. *Am. J. Physiol. Gastrointest. Liver Physiol.* 279: G775-G780.
- Murakami, Y., Kohyama, N., Kobayashi, Y., Ohbayashi, M., Ohtani, H., Sawada, Y. and Yamamoto, T. 2005. Functional characterization of human monocarboxylate transporter 6 (slc16a5). *Drug Metab. Dispos.* 33: 1845-1851.
- Bonen, A., Heynen, M. and Hatta, H. 2006. Distribution of monocarboxylate transporters MCT1-MCT8 in rat tissues and human skeletal muscle. *Appl. Physiol. Nutr. Metab.* 31: 31-39.
- Visser, W.E., Friesema, E.C., Jansen, J. and Visser, T.J. 2007. Thyroid hormone transport by monocarboxylate transporters. *Best Pract. Res. Clin. Endocrinol. Metab.* 21: 223-236.

CHROMOSOMAL LOCATION

Genetic locus: SLC16A5 (human) mapping to 17q25.1.

PRODUCT

MCT6 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see MCT6 shRNA Plasmid (h): sc-106212-SH and MCT6 shRNA (h) Lentiviral Particles: sc-106212-V as alternate gene silencing products.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

MCT6 siRNA (h) is recommended for the inhibition of MCT6 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor MCT6 gene expression knockdown using RT-PCR Primer: MCT6 (h)-PR: sc-106212-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.