

mucolipin 3 siRNA (h): sc-106264

BACKGROUND

The mucolipins are a family of cation channel proteins designated mucolipin 1, mucolipin 2 and mucolipin 3. Mucolipins are predominantly expressed within the endocytic pathway and may regulate membrane traffic and/or degradation. Mucolipin 3, also known as TRPML3, TRP-ML3 or MCOLN3, is a 553 amino acid multi-pass membrane protein that belongs to the Polycystin family of transient receptors and localizes to lysosomes. Also expressed in hair cells, mucolipin 3 may play a critical role at the ankle-link region during hair-bundle growth. Mutations in mucolipin 3 have been linked to hair cell degeneration leading to hearing loss and balance disorders. Defects in mucolipin 3 may be associated with mucopolipidosis type IV, an autosomal recessive lysosomal storage disorder that causes severe neurodegeneration, achlorhydria and visual impairments. Two isoforms exist due to alternative splicing events.

REFERENCES

1. Song, Y., et al. 2006. TRPML cation channels regulate the specialized lysosomal compartment of vertebrate B-lymphocytes. *Eur. J. Cell Biol.* 85: 1253-1264.
2. Venkatachalam, K., et al. 2006. Lysosomal localization of TRPML3 depends on TRPML2 and the mucopolipidosis-associated protein TRPML1. *J. Biol. Chem.* 281: 17517-17527.
3. Cuajungco, M.P., et al. 2007. TRP channels as candidates for hearing and balance abnormalities in vertebrates. *Biochim. Biophys. Acta* 1772: 1022-1027.
4. Atiba-Davies, M., et al. 2007. TRPML3 and hearing loss in the varitint-waddler mouse. *Biochim. Biophys. Acta* 1772: 1028-1031.
5. Kim, H.J., et al. 2007. Gain-of-function mutation in TRPML3 causes the mouse varitint-waddler phenotype. *J. Biol. Chem.* 282: 36138-36142.
6. Karacsonyi, C., et al. 2007. Mucolipin-2 localizes to the Arf6-associated pathway and regulates recycling of GPI-APs. *Traffic* 8: 1404-1414.
7. van Aken, A.F., et al. 2008. TRPML3 mutations cause impaired mechano-electrical transduction and depolarization by an inward-rectifier cation current in auditory hair cells of varitint-waddler mice. *J. Physiol.* 586: 5403-5418.

CHROMOSOMAL LOCATION

Genetic locus: MCOLN3 (human) mapping to 1p22.3.

PRODUCT

mucolipin 3 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see mucolipin 3 shRNA Plasmid (h): sc-106264-SH and mucolipin 3 shRNA (h) Lentiviral Particles: sc-106264-V as alternate gene silencing products.

For independent verification of mucolipin 3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-106264A, sc-106264B and sc-106264C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

mucolipin 3 siRNA (h) is recommended for the inhibition of mucolipin 3 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contains a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor mucolipin 3 gene expression knockdown using RT-PCR Primer: mucolipin 3 (h)-PR: sc-106264-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.