



NBPF1 siRNA (h): sc-106281

BACKGROUND

NBPF1 (neuroblastoma breakpoint family, member 1), also known as AD2, NBG, AB13, AB14 or AB23, is a 1,214 amino acid protein that localizes to the cytoplasm and contains eight NBPF domains. Expressed in a variety of tissues, NBPF1 belongs to the neuroblastoma breakpoint family and is encoded by a gene which maps to human chromosome 1. Chromosome 1 spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome. Chromosome 1 houses a large number of disease-associated genes, including those that are involved in familial adenomatous polyposis, Stickler syndrome, Parkinson's disease, Gaucher disease, schizophrenia and Usher syndrome. Aberrations in chromosome 1 are found in a variety of cancers, including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

1. Weise, A., Starke, H., Mrasek, K., Claussen, U. and Liehr, T. 2005. New insights into the evolution of chromosome 1. *Cytogenet. Genome Res.* 108: 217-222.
2. Vandepoele, K., Van Roy, N., Staes, K., Speleman, F. and van Roy, F. 2005. A novel gene family NBPF: intricate structure generated by gene duplications during primate evolution. *Mol. Biol. Evol.* 22: 2265-2274.
3. Marzin, Y., Jamet, D., Douet-Guilbert, N., Morel, F., Le Bris, M.J., Morice, P., Abgrall, J.F., Berthou, C. and De Braekeleer, M. 2006. Chromosome 1 abnormalities in multiple myeloma. *Anticancer Res.* 26: 953-959.
4. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 610501. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Vandepoele, K., Andries, V., Van Roy, N., Staes, K., Vandesompele, J., Laureys, G., De Smet, E., Berx, G., Speleman, F. and van Roy, F. 2008. A constitutional translocation t(1;17)(p36.2;q11.2) in a neuroblastoma patient disrupts the human NBPF1 and ACCN1 genes. *PLoS ONE* 3: e2207.
6. Vandepoele, K., Andries, V. and van Roy, F. 2009. The NBPF1 promoter has been recruited from the unrelated EVI5 gene before simian radiation. *Mol. Biol. Evol.* 26: 1321-1332.

CHROMOSOMAL LOCATION

Genetic locus: NBPF1 (human) mapping to 1p36.13.

PRODUCT

NBPF1 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see NBPF1 shRNA Plasmid (h): sc-106281-SH and NBPF1 shRNA (h) Lentiviral Particles: sc-106281-V as alternate gene silencing products.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

NBPF1 siRNA (h) is recommended for the inhibition of NBPF1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor NBPF1 gene expression knockdown using RT-PCR Primer: NBPF1 (h)-PR: sc-106281-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.