

NGL-1 siRNA (h): sc-106302

BACKGROUND

NGL-1 (netrin G1 ligand), also known as leucine-rich repeat-containing protein 4C, is a single pass type I membrane protein that acts as a cell adhesion molecule. It contains nine leucine-rich repeats (LRR) and one Ig-like C2-type domain. NGL-1 is predominantly expressed in the striatum and the cerebral cortex of both the embryonic and adult brain. NGL-1 specifically interacts with netrin G1 (a molecule involved in axon guidance in the developing central nervous system) via its LRR region. NGL-1 plays a role in the regulation of neurite outgrowth of developing thalamic neurons. Soluble NGL-1 inhibits thalamic axon outgrowth while NGL-1 that is bound to the surface of developing thalamocortical axons stimulates growth. NGL-1 also interacts with Whirlin possibly stabilizing interstereociliar links.

REFERENCES

1. Lin, J.C., Ho, W.H., Gurney, A. and Rosenthal, A. 2003. The netrin-G1 ligand NGL-1 promotes the outgrowth of thalamocortical axons. *Nat. Neurosci.* 6: 1270-1276.
2. Delprat, B., Michel, V., Goodyear, R., Yamasaki, Y., Michalski, N., El-Amraoui, A., Perfettini, I., Legrain, P., Richardson, G., Hardelin, J.P. and Petit, C. 2005. Myosin XVa and whirlin, two deafness gene products required for hair bundle growth, are located at the stereocilia tips and interact directly. *Hum. Mol. Genet.* 14: 401-410.
3. Miyashita, T., Nishimura-Akiyoshi, S., Itohara, S. and Rockland, K.S. 2005. Strong expression of netrin G2 in the monkey claustrum. *Neuroscience* 136: 487-496.
4. Meerabux, J.M., Ohba, H., Fukasawa, M., Suto, Y., Aoki-Suzuki, M., Nakashiba, T., Nishimura, S., Itohara, S. and Yoshikawa, T. 2005. Human netrin G1 isoforms show evidence of differential expression. *Genomics* 86: 112-116.
5. Morimura, N., Inoue, T., Katayama, K. and Aruga, J. 2006. Comparative analysis of structure, expression and PSD95-binding capacity of Lfn, a novel family of neuronal transmembrane proteins. *Gene* 380: 72-83.
6. Ko, J., Kim, S., Chung, H.S., Kim, K., Han, K., Kim, H., Jun, H., Kaang, B.K. and Kim, E. 2006. SALM synaptic cell adhesion-like molecules regulate the differentiation of excitatory synapses. *Neuron* 50: 233-245.
7. Amos, C.I., Chen, W.V., Lee, A., Li, W., Kern, M., Lundsten, R., Batliwalla, F., Wener, M., Remmers, E., Kastner, D.A., Criswell, L.A., Seldin, M.F. and Gregersen, P.K. 2006. High-density SNP analysis of 642 Caucasian families with rheumatoid arthritis identifies two new linkage regions on 11p12 and 2q33. *Genes Immun.* 7: 277-286.
8. Chen, Y., Aulia, S., Li, L. and Tang, B.L. 2006. AMIGO and friends: an emerging family of brain-enriched, neuronal growth modulating, type I transmembrane proteins with leucine-rich repeats (LRR) and cell adhesion molecule motifs. *Brain Res. Rev.* 51: 265-274.
9. Nishimura-Akiyoshi, S., Niimi, K., Nakashiba, T. and Itohara, S. 2007. Axonal netrin-Gs transneurally determine lamina-specific subdendritic segments. *Proc. Natl. Acad. Sci. USA* 104: 14801-14806.

CHROMOSOMAL LOCATION

Genetic locus: LRRC4C (human) mapping to 11p12.

PRODUCT

NGL-1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see NGL-1 shRNA Plasmid (h): sc-106302-SH and NGL-1 shRNA (h) Lentiviral Particles: sc-106302-V as alternate gene silencing products.

For independent verification of NGL-1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-106302A, sc-106302B and sc-106302C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

NGL-1 siRNA (h) is recommended for the inhibition of NGL-1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor NGL-1 gene expression knockdown using RT-PCR Primer: NGL-1 (h)-PR: sc-106302-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.