

Nkx-2.4 siRNA (h): sc-106305

BACKGROUND

Members of the NK-2 family of homeodomain proteins, which include Nkx-2.2, Nkx-2.4, Nkx-2.5, Nkx-2.6 and many others, are key regulators of growth and development in several tissues, including brain, heart and pancreas. Nkx-2.2 is responsible for directing ventral neuronal patterning in response to graded Shh signaling. Nkx-2.4 is a 354 amino acid nuclear protein that is detected in the posterior hypothalamus of the embryo and also in adult testis. Nkx-2.5, also designated cardiac specific homeobox protein (Csx), is a homolog of the *Drosophila* tinman protein and is essential for normal cardiovascular development. Nkx-2.6, also a homolog of the *Drosophila* tinman protein, is expressed in the caudal pharyngeal pouches, the caudal heart progenitors, the sinus venosus, the outflow tract of the heart and in a short segment of the gut between stages E8.5 and E10.5 of embryogenesis.

REFERENCES

1. Price, M., Lazzaro, D., Pohl, T., Mattei, M.G., Rüther, U., Olivo, J.C., Duboule, D. and Di Lauro, R. 1992. Regional expression of the homeobox gene Nkx-2.2 in the developing mammalian forebrain. *Neuron* 8: 241-255.
2. Sussel, L., Kalamaras, J., Hartigan-O'Connor, D.J., Meneses, J.J., Pedersen, R.A., Rubenstein, J.L. and German, M.S. 1998. Mice lacking the homeodomain transcription factor Nkx2.2 have diabetes due to arrested differentiation of pancreatic beta cells. *Development* 125: 2213-2221.
3. Pabst, O., Herbrand, H. and Arnold, H.H. 1998. Nkx2-9 is a novel homeobox transcription factor which demarcates ventral domains in the developing mouse CNS. *Mech. Dev.* 73: 85-93.
4. McMahon, A.P. 2000. Neural patterning: the role of Nkx genes in the ventral spinal cord. *Genes Dev.* 14: 2261-2264.
5. Wang, C.C., Brodnicki, T., Copeland, N.G., Jenkins, N.A. and Harvey, R.P. 2000. Conserved linkage of NK-2 homeobox gene pairs Nkx2-2/2-4 and Nkx2-1/2-9 in mammals. *Mamm. Genome* 11: 466-468.
6. Holland, P.W., Booth, H.A. and Bruford, E.A. 2007. Classification and nomenclature of all human homeobox genes. *BMC Biol.* 5: 47.
7. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 607808. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: NKX2-4 (human) mapping to 20p11.22.

PRODUCT

Nkx-2.4 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Nkx-2.4 shRNA Plasmid (h): sc-106305-SH and Nkx-2.4 shRNA (h) Lentiviral Particles: sc-106305-V as alternate gene silencing products.

For independent verification of Nkx-2.4 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-106305A, sc-106305B and sc-106305C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Nkx-2.4 siRNA (h) is recommended for the inhibition of Nkx-2.4 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Nkx-2.4 gene expression knockdown using RT-PCR Primer: Nkx-2.4 (h)-PR: sc-106305-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.