

SSBP3 siRNA (m): sc-106565

BACKGROUND

The single-stranded-DNA-binding proteins (SSBs) are essential for DNA function in prokaryotic and eukaryotic cells, as well as in mitochondria, bacteria and viruses. SSBP3 (single stranded DNA binding protein 3), also known as CSDP or SSDP, is a 388 amino acid protein that localizes to the nucleus and contains one LisH domain. Highly expressed in spleen, bone marrow, thymus, lymph node, kidney, brain, heart and skeletal muscle, SSBP3 binds to the single-stranded polypyrimidine sequences in the promotor of COLA2 genes and is thought to regulate the transcription of COLA2 genes. SSBP3 is expressed as multiple alternatively spliced isoforms that are encoded by a gene which maps to human chromosome 1.

REFERENCES

1. Raval-Fernandes, S., Kickhoefer, V.A. and Rome, L.H. 1999. Cloning of a cDNA encoding a sequence-specific single-stranded-DNA-binding protein from *Rattus norvegicus*. *Gene* 237: 201-207.
2. Castro, P., Liang, H., Liang, J.C. and Nagarajan, L. 2002. A novel, evolutionarily conserved gene family with putative sequence-specific single-stranded DNA-binding activity. *Genomics* 80: 78-85.
3. Chen, L., Segal, D., Hukriede, N.A., Podtelejnikov, A.V., Bayarsaihan, D., Kennison, J.A., Ogryzko, V.V., Dawid, I.B. and Westphal, H. 2002. Ssdp proteins interact with the LIM-domain-binding protein Ldb1 to regulate development. *Proc. Natl. Acad. Sci. USA* 99: 14320-14325.
4. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607390. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. van Meyel, D.J., Thomas, J.B. and Agulnick, A.D. 2003. Ssdp proteins bind to LIM-interacting co-factors and regulate the activity of LIM-homeodomain protein complexes *in vivo*. *Development* 130: 1915-1925.
6. Wu, L. 2006. Structure and functional characterization of single-strand DNA binding protein SSDP1: carboxyl-terminal of SSDP1 has transcription activity. *Biochem. Biophys. Res. Commun.* 339: 977-984.
7. Xu, Z., Meng, X., Cai, Y., Liang, H., Nagarajan, L. and Brandt, S.J. 2007. Single-stranded DNA-binding proteins regulate the abundance of LIM domain and LIM domain-binding proteins. *Genes Dev.* 21: 942-955.
8. Cai, Y., Xu, Z., Nagarajan, L. and Brandt, S.J. 2008. Single-stranded DNA-binding proteins regulate the abundance and function of the LIM-homeodomain transcription factor LHX2 in pituitary cells. *Biochem. Biophys. Res. Commun.* 373: 303-308.

CHROMOSOMAL LOCATION

Genetic locus: Ssbp3 (mouse) mapping to 4 C7.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

SSBP3 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SSBP3 shRNA Plasmid (m): sc-106565-SH and SSBP3 shRNA (m) Lentiviral Particles: sc-106565-V as alternate gene silencing products.

For independent verification of SSBP3 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-106565A, sc-106565B and sc-106565C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SSBP3 siRNA (m) is recommended for the inhibition of SSBP3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SSBP3 gene expression knockdown using RT-PCR Primer: SSBP3 (m)-PR: sc-106565-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.