

TAF II p20 siRNA (m): sc-106595

BACKGROUND

TFIID is a general transcription factor that initiates preinitiation complex assembly through direct interaction with the TATA promoter element. Functioning as a multisubunit complex consisting of a small TATA-binding polypeptide and other TBP-associated factors (TAFs), TFIID mediates promoter responses to various transcriptional activators and repressors. TAF II p20, also known as TAF12, TAF15, TAF2J or TAFII20, is a 161 amino acid subunit of TFIID that localizes to the nucleus and contains one histone-fold domain. Expressed ubiquitously, TAF II p20 interacts with other members of the TFIID complex and, via this interaction, plays a role in mediating transcriptional activation and repression. Two isoforms of TAF II p20 exist due to alternative splicing events.

REFERENCES

1. Holstege, F.C., Fiedler, U. and Timmers, H.T. 1997. Three transitions in the RNA polymerase II transcription complex during initiation. *EMBO J.* 16: 7468-7480.
2. Ogryzko, V.V., Kotani, T., Zhang, X., Schiltz, R.L., Howard, T., Yang, X.J., Howard, B.H., Qin, J. and Nakatani, Y. 1998. Histone-like TAFs within the PCAF histone acetylase complex. *Cell* 94: 35-44.
3. Gangloff, Y.G., Werten, S., Romier, C., Carre, L., Poch, O., Moras, D. and Davidson, I. 2000. The human TFIID components TAF_{II}135 and TAF_{II}20 and the yeast SAGA components ADA1 and TAF_{II}68 heterodimerize to form histone-like pairs. *Mol. Cell. Biol.* 20: 340-351.
4. Werten, S., Mitschler, A., Romier, C., Gangloff, Y.G., Thuault, S., Davidson, I. and Moras, D. 2002. Crystal structure of a subcomplex of human transcription factor TFIID formed by TATA binding protein-associated factors hTAF4 (hTAF_{II}135) and hTAF12 (hTAF_{II}20). *J. Biol. Chem.* 277: 45502-45509.
5. Guermah, M., Ge, K., Chiang, C.M. and Roeder, R.G. 2003. The TBN protein, which is essential for early embryonic mouse development, is an inducible TAFII implicated in adipogenesis. *Mol. Cell* 12: 991-1001.
6. Cavusoglu, N., Brand, M., Tora, L. and Van Dorsselaer, A. 2003. Novel subunits of the TATA binding protein free TAFII-containing transcription complex identified by matrix-assisted laser desorption/ionization-time of flight mass spectrometry following one-dimensional gel electrophoresis. *Proteomics* 3: 217-223.

CHROMOSOMAL LOCATION

Genetic locus: TAF12 (mouse) mapping to 4 D2.3.

PRODUCT

TAF II p20 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TAF II p20 shRNA Plasmid (m): sc-106595-SH and TAF II p20 shRNA (m) Lentiviral Particles: sc-106595-V as alternate gene silencing products.

For independent verification of TAF II p20 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-106595A, sc-106595B and sc-106595C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TAF II p20 siRNA (m) is recommended for the inhibition of TAF II p20 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

TAF II p20 (B-6): sc-514619 is recommended as a control antibody for monitoring of TAF II p20 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TAF II p20 gene expression knockdown using RT-PCR Primer: TAF II p20 (m)-PR: sc-106595-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.