

Sec61 β siRNA (h): sc-106868

BACKGROUND

The Sec61 complex forms the core element of the protein translation complex in the rough endoplasmic reticulum membrane. The complex also associates with two ubiquitous ER membrane proteins: Sec62 (also designated human translocation protein 1 or HTP1) and Sec63. The complex forms a two-way channel that transports proteins both into the ER and back to the cytosol for degradation. Specifically, it appears the β subunit facilitates the escort of proteins back to the cytoplasm for degradation by the proteasome or by other proteolytic systems.

REFERENCES

- Rapoport, T.A., et al. 1996. Protein transport across the eukaryotic endoplasmic reticulum and bacterial inner membranes. *Annu. Rev. Biochem.* 65: 271-303.
- Daimon, M., et al. 1997. Identification of a human cDNA homologue to the *Drosophila* translocation protein 1 (Dtrp1). *Biochem. Biophys. Res. Commun.* 230: 100-104.
- Bebok, Z., et al. 1998. The mechanism underlying cystic fibrosis transmembrane conductance regulator transport from the endoplasmic reticulum to the proteasome includes Sec61 β and a cytosolic, deglycosylated intermediary. *J. Biol. Chem.* 273: 29873-29878.
- Romisch, K. 1999. Surfing the Sec61 channel: bidirectional protein translocation across the ER membrane. *J. Cell Sci.* 112: 4185-4191.
- Raden, D., et al. 2000. Role of the cytoplasmic segments of Sec61 α in the ribosome-binding and translocation-promoting activities of the Sec61 complex. *J. Cell Biol.* 150: 53-64.
- Meyer, H.A., et al. 2000. Mammalian Sec61 is associated with Sec62 and Sec63. *J. Biol. Chem.* 275: 14550-14557.
- Levy, R., et al. 2001. *In vitro* binding of ribosomes to the β subunit of the Sec61p protein translocation complex. *J. Biol. Chem.* 276: 2340-2346.
- Alder, N.N., et al. 2005. The molecular mechanisms underlying BiP-mediated gating of the Sec61 translocon of the endoplasmic reticulum. *J. Cell Biol.* 168: 389-399.
- Wilson, C.M., et al. 2005. Ribophorin I associates with a subset of membrane proteins after their integration at the sec61 translocon. *J. Biol. Chem.* 280: 4195-4206.

CHROMOSOMAL LOCATION

Genetic locus: SEC61B (human) mapping to 9q22.33.

PRODUCT

Sec61 β siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Sec61 β shRNA Plasmid (h): sc-106868-SH and Sec61 β shRNA (h) Lentiviral Particles: sc-106868-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20 $^{\circ}$ C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20 $^{\circ}$ C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Sec61 β siRNA (h) is recommended for the inhibition of Sec61 β expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Sec61 β (E-6): sc-393633 is recommended as a control antibody for monitoring of Sec61 β gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz MarkerTM Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.