



TNP1 siRNA (h): sc-106912

BACKGROUND

TNP1 (transition protein 1), also known as TP1 or STP1 (spermatid nuclear transition protein 1), is a 55 amino acid nuclear protein that belongs to the nuclear transition protein family. Expressed specifically in the testis, TNP1 is a spermatid-specific protein that functions as a basic chromosomal transition protein during the conversion of nucleosomal chromatin to the compact, non-nucleosomal form found in sperm nuclei. During spermatogenesis, TNP1, a product of the haploid genome, replaces histone and is subsequently replaced by protamines in mature sperm. Although TNP1 is present for only a short amount of time during sperm maturation, it is thought to play an essential role in male fertility.

REFERENCES

1. Luerssen, H., Mattei, M.G., Schröter, M., Grzeschik, K.H., Adham, I.M. and Engel, W. 1990. Nucleotide sequence of the gene for human transition protein 1 and its chromosomal localization on chromosome 2. *Genomics* 8: 324-330.
2. Yelick, P.C., Kozak, C., Kwon, Y.K., Seldin, M.F. and Hecht, N.B. 1991. The mouse transition protein 1 gene contains a B1 repetitive element and is located on chromosome 1. *Genomics* 11: 687-694.
3. Yu, Y.E., Zhang, Y., Unni, E., Shirley, C.R., Deng, J.M., Russell, L.D., Weil, M.M., Behringer, R.R. and Meistrich, M.L. 2000. Abnormal spermatogenesis and reduced fertility in transition nuclear protein 1-deficient mice. *Proc. Natl. Acad. Sci. USA* 97: 4683-4688.
4. Meistrich, M.L., Mohapatra, B., Shirley, C.R. and Zhao, M. 2003. Roles of transition nuclear proteins in spermiogenesis. *Chromosoma* 111: 483-488.
5. Miyagawa, Y., Nishimura, H., Tsujimura, A., Matsuoka, Y., Matsumiya, K., Okuyama, A., Nishimune, Y. and Tanaka, H. 2005. Single-nucleotide polymorphisms and mutation analyses of the TNP1 and TNP2 genes of fertile and infertile human male populations. *J. Androl.* 26: 779-786.
6. Tüttelmann, F., Rajpert-De Meyts, E., Nieschlag, E. and Simoni, M. 2007. Gene polymorphisms and male infertility—a meta-analysis and literature review. *Reprod. Biomed. Online* 15: 643-658.
7. Okada, Y., Scott, G., Ray, M.K., Mishina, Y. and Zhang, Y. 2007. Histone demethylase JHDM2A is critical for TNP1 and Prm1 transcription and spermatogenesis. *Nature* 450: 119-123.
8. Mitchell, V., Lefebvre-Khalil, V., Thomas, P., Rigot, J.M. and Steger, K. 2007. Transition protein 1 mRNA expression is not related to pregnancy rate in azoospermic men undergoing TESE-ICSI. *Andrologia* 39: 124-127.
9. Jedrzejczak, P., Kempisty, B., Bryja, A., Mostowska, M., Depa-Martynow, M., Pawelczyk, L. and Jagodzinski, P.P. 2007. Quantitative assessment of transition proteins 1, 2 spermatid-specific linker Histone H1-like protein transcripts in spermatozoa from normozoospermic and asthenozoospermic men. *Arch. Androl.* 53: 199-205.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: TNP1 (human) mapping to 2q35.

PRODUCT

TNP1 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TNP1 shRNA Plasmid (h): sc-106912-SH and TNP1 shRNA (h) Lentiviral Particles: sc-106912-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TNP1 siRNA (h) is recommended for the inhibition of TNP1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.