

# POMZP3 siRNA (h): sc-106945

## BACKGROUND

POMZP3 (POM and ZP3 fusion protein) is a 210 amino acid fusion protein expressed in a variety of human tissues, including colon, spleen, testis, thymus, ovary, lymphocytes, pancreas and small intestine. POMZP3 is encoded by a gene that resulted from the fusion of duplicated POM121 and ZP3 DNA sequences. The 5' end of POMZP3 shares 77% identity with POM121, while the 3' end shares 99% identity with ZP3. POM121 is an integral nuclear pore membrane protein that may contribute to the permeability of the nuclear pore complex (NPC). Unlike POM121, POMZP3 does not contain the nuclear pore localization motif and, therefore, most likely does not associate with the NPC. ZP3 functions as a sperm receptor ligand of the zona pellucida and plays a key role in the first step of the fertilization process. Following sperm-oocyte binding, ZP3 triggers the sperm acrosome reaction that releases the protein machinery, enabling a spermatozoon to penetrate the zona pellucida.

## REFERENCES

1. Liang, L.F., et al. 1990. Oocyte-specific expression of mouse Zp-2: developmental regulation of the zona pellucida genes. *Mol. Cell. Biol.* 10: 1507-1515.
2. Dean, J. 1992. Biology of mammalian fertilization: role of the zona pellucida. *J. Clin. Invest.* 89: 1055-1059.
3. Kipersztok, S., et al. 1995. POM-ZP3, a bipartite transcript derived from human ZP3 and POM121 homologue. *Genomics* 25: 354-359.
4. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 600587. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Gupta, S.K., et al. 2003. Zona pellucida glycoproteins based immunocontraceptive vaccines: strategies for development and their applications. *Indian J. Exp. Biol.* 41: 682-693.
6. Jazwinska, A. and Affolter, M. 2004. A family of genes encoding zona pellucida (ZP) domain proteins is expressed in various epithelial tissues during *Drosophila* embryogenesis. *Gene Expr. Patterns* 4: 413-421.
7. Wassarman, P.M., et al. 2004. Mouse zona pellucida genes and glycoproteins. *Cytogenet. Genome Res.* 105: 228-234.
8. Wassarman, P.M., et al. 2004. Egg-sperm interactions at fertilization in mammals. *Eur. J. Obstet. Gynecol. Reprod. Biol.* 115: S57-S60.
9. Naz, R.K. and Rajesh, C. 2005. Gene knockouts that cause female infertility: search for novel contraceptive targets. *Front. Biosci.* 10: 2447-2459.

## CHROMOSOMAL LOCATION

Genetic locus: POMZP3 (human) mapping to 7q11.23.

## PRODUCT

POMZP3 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see POMZP3 shRNA Plasmid (h): sc-106945-SH and POMZP3 shRNA (h) Lentiviral Particles: sc-106945-V as alternate gene silencing products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

POMZP3 siRNA (h) is recommended for the inhibition of POMZP3 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

POMZP3 (YY-47): sc-100902 is recommended as a control antibody for monitoring of POMZP3 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor POMZP3 gene expression knockdown using RT-PCR Primer: POMZP3 (h)-PR: sc-106945-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.