

3 β -HSD2 siRNA (m): sc-108885

BACKGROUND

3 β -hydroxysteroid dehydrogenase (3 β -HSD), also known as HSD3B1 or HSD3B3, is a bifunctional enzyme that plays a crucial role in the synthesis of all classes of hormonal steroids. Two human 3 β -HSD proteins, designated type I (3 β -HSD) and type II (3 β -HSD2), are expressed by different genes and function in different areas of the body. Localized to the membrane of the endoplasmic reticulum (ER) and expressed in testis, ovaries and adrenal gland, 3 β -HSD2 is the type II protein that catalyzes the oxidative conversion of Δ^5 -ene-3 β -hydroxysteroid, as well as the conversion of various ketosteroids. Defects in the gene encoding 3 β -HSD2 are the cause of adrenal hyperplasia type 2 (AH2), a form of recessive congenital adrenal hyperplasia that is characterized by excess androgen which can lead to ambiguous genitalia and rapid somatic growth.

REFERENCES

1. Thomas, J.L., et al. 2002. Structure/function relationships responsible for the kinetic differences between human type 1 and type 2 3 β -hydroxysteroid dehydrogenase and for the catalysis of the type 1 activity. *J. Biol. Chem.* 277: 42795-42801.
2. Thomas, J.L., et al. 2003. Structure/function relationships responsible for coenzyme specificity and the isomerase activity of human type 1 3 β -hydroxysteroid dehydrogenase/isomerase. *J. Biol. Chem.* 278: 35483-35490.
3. Foti, D.M., et al. 2004. YY1 binding within the human HSD3B2 gene intron 1 is required for maximal basal promoter activity: identification of YY1 as the 3 β 1-A factor. *J. Mol. Endocrinol.* 33: 99-119.
4. Thomas, J.L., et al. 2004. Serine 124 completes the Tyr, Lys and Ser triad responsible for the catalysis of human type 1 3 β -hydroxysteroid dehydrogenase. *J. Mol. Endocrinol.* 33: 253-261.
5. Carbanaru, G., et al. 2004. The hormonal phenotype of nonclassic 3 β -hydroxysteroid dehydrogenase (HSD3B) deficiency in hyperandrogenic females is associated with Insulin-resistant polycystic ovary syndrome and is not a variant of inherited HSD3B2 deficiency. *J. Clin. Endocrinol. Metab.* 89: 783-794.
6. Thomas, J.L., et al. 2007. Structure/function of human type 1 3 β -hydroxysteroid dehydrogenase: an intrasubunit disulfide bond in the Rossmann-fold domain and a Cys residue in the active site are critical for substrate and coenzyme utilization. *J. Steroid Biochem. Mol. Biol.* 107: 80-87.
7. Wang, L., et al. 2007. Human 3 β -hydroxysteroid dehydrogenase types 1 and 2: gene sequence variation and functional genomics. *J. Steroid Biochem. Mol. Biol.* 107: 88-99.
8. Park, J.Y., et al. 2007. Association between polymorphisms in HSD3B1 and UGT2B17 and prostate cancer risk. *Urology* 70: 374-379.
9. Xing, Y., et al. 2008. The farnesoid X receptor regulates transcription of 3 β -hydroxysteroid dehydrogenase type 2 in human adrenal cells. *Mol. Cell. Endocrinol.* 299: 153-162.

CHROMOSOMAL LOCATION

Genetic locus: Hsd3b2 (mouse) mapping to 3 F2.2.

PRODUCT

3 β -HSD2 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see 3 β -HSD2 shRNA Plasmid (m): sc-108885-SH and 3 β -HSD2 shRNA (m) Lentiviral Particles: sc-108885-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

3 β -HSD2 siRNA (m) is recommended for the inhibition of 3 β -HSD2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor 3 β -HSD2 gene expression knockdown using RT-PCR Primer: 3 β -HSD2 (m)-PR: sc-108885-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.