

HspBP1 (h): 293 Lysate: sc-110611

BACKGROUND

HSP 70-interacting protein (HspBP1) belongs to a family of eukaryotic proteins identified as nucleotide exchange factors for HSP 70, which exhibit varying degrees of compartment and species specificity. HspBP1 interferes with the CHIP-induced degradation of immature forms of the cystic fibrosis transmembrane conductance regulator (CFTR) and stimulates CFTR maturation. HspBP1 binds to HSP 70, inhibits its activity and promotes dissociation of nucleotides from the HSP 70 ATPase domain. It is a protein mainly expressed in heart and skeletal muscle.

REFERENCES

1. Raynes, D.A. and Guerriero, V. 2000. Isolation and characterization of isoforms of HspBP1, inhibitors of HSP 70. *Biochim. Biophys. Acta* 1490: 203-207.
2. Kabani, M., McLellan, C., Raynes, D.A., Guerriero, V. and Brodsky, J.L. 2002. HspBP1, a homologue of the yeast Fes1 and Sls1 proteins, is an HSC 70 nucleotide exchange factor. *FEBS Lett.* 531: 339-342.
3. Raynes, D.A., Graner, M.W., Bagatell, R., McLellan, C. and Guerriero, V. 2003. Increased expression of the HSP 70 cochaperone HspBP1 in tumors. *Tumour Biol.* 24: 281-285.
4. McLellan, C.A., Raynes, D.A. and Guerriero, V. 2003. HspBP1, an HSP 70 cochaperone, has two structural domains and is capable of altering the conformation of the HSP 70 ATPase domain. *J. Biol. Chem.* 278: 19017-19022.
5. Tanimura, S. and Kohno, M. 2004. Heat shock protein 70 binding protein 1 induces enhanced apoptotic response against anticancer drugs in tumor cells. *Nippon Rinsho* 62: 1291-1296.
6. Alberti, S., Böhse, K., Arndt, V., Schmitz, A. and Höhfeld, J. 2004. The cochaperone HspBP1 inhibits the CHIP ubiquitin ligase and stimulates the maturation of the cystic fibrosis transmembrane conductance regulator. *Mol. Biol. Cell* 15: 4003-4010.
7. Papp, D., Prohászka, Z., Kocsis, J., Füst, G., Bánhegyi, D., Raynes, D.A. and Guerriero, V. 2005. Development of a sensitive assay for the measurement of antibodies against heat shock protein binding protein 1 (HspBP1): increased levels of anti-HspBP1 IgG are prevalent in HIV infected subjects. *J. Med. Virol.* 76: 464-469.
8. Shomura, Y., Dragovic, Z., Chang, H.C., Tzvetkov, N., Young, J.C., Brodsky, J.L., Guerriero, V., Hartl, F.U. and Bracher, A. 2005. Regulation of HSP 70 function by HspBP1: structural analysis reveals an alternate mechanism for HSP 70 nucleotide exchange. *Mol. Cell* 17: 367-379.

CHROMOSOMAL LOCATION

Genetic locus: HSPBP1 (human) mapping to 19q13.42.

PRODUCT

HspBP1 (h): 293 Lysate represents a lysate of human HspBP1 transfected 293 cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

HspBP1 (h): 293 Lysate is suitable as a Western Blotting positive control for human reactive HspBP1 antibodies. Recommended use: 10-20 µl per lane.


Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

HspBP1 (D-2): sc-166315 is recommended as a positive control antibody for Western Blot analysis of enhanced human HspBP1 expression in HspBP1 transfected 293 cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



HspBP1 (D-2): sc-166315. Western blot analysis of HspBP1 expression in non-transfected: sc-110760 (A) and human HspBP1 transfected: sc-110611 (B) 293 whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.