

Galactose Mutarotase (h): 293 Lysate: sc-112214

BACKGROUND

Galactose Mutarotase is a member of the aldose epimerase family and is involved in hexose metabolism. Through its catalytic activity, Galactose Mutarotase converts β -aldose to α -aldose on several sugars, including D-glucose, L-arabinose and D-xylose. Found in the cytoplasm of most cells, Galactose Mutarotase plays a key role in galactose metabolism by catalyzing the conversion of β -D-galactose to α -D-galactose. The enzyme contains two residues, Glu 304 and His 170, that are critical for catalysis, as well as His 96 and Asp 243, which are important for proper substrate recognition by the active site. No known diseases have been associated with mutations in the Galactose Mutarotase gene, although inhibition of Galactose Mutarotase activity could potentially be associated with a buildup of unmetabolized sugars during metabolism.

REFERENCES

1. Beebe, J.A. and Frey, P.A. 1998. Galactose Mutarotase: purification, characterization, and investigations of two important histidine residues. *Biochemistry* 37: 14989-14997.
2. Beebe, J.A., Arabshahi, A., Clifton, J.G., Ringe, D., Petsko, G.A. and Frey, P.A. 2003. Galactose Mutarotase: pH dependence of enzymatic mutarotation. *Biochemistry* 42: 4414-4420.
3. Thoden, J.B., Kim, J., Raushel, F.M. and Holden, H.M. 2003. The catalytic mechanism of Galactose Mutarotase. *Protein Sci.* 12: 1051-1059.
4. Thoden, J.B., Timson, D.J., Reece, R.J. and Holden, H.M. 2004. Molecular structure of human Galactose Mutarotase. *J. Biol. Chem.* 279: 23431-23437.
5. Kim, I., Kim, E., Yoo, S., Shin, D., Min, B., Song, J. and Park, C. 2004. Ribose utilization with an excess of mutarotase causes cell death due to accumulation of methylglyoxal. *J. Bacteriol.* 186: 7229-7235.
6. Thoden, J.B. and Holden, H.M. 2005. The molecular architecture of Galactose Mutarotase/UDP-galactose 4-epimerase from *Saccharomyces cerevisiae*. *J. Biol. Chem.* 280: 21900-21907.
7. Ryu, K.S., Kim, J.I., Cho, S.J., Park, D., Park, C., Cheong, H.K., Lee, J.O. and Choi, B.S. 2005. Structural insights into the monosaccharide specificity of *Escherichia coli* rhamnose mutarotase. *J. Mol. Biol.* 349: 153-162.
8. Barreto, M., Jedlicki, E. and Holmes, D.S. 2005. Identification of a gene cluster for the formation of extracellular polysaccharide precursors in the chemolithoautotroph *Acidithiobacillus ferrooxidans*. *Appl. Environ. Microbiol.* 71: 2902-2909.
9. Miron, S., Réfregiers, M., Gilles, A.M. and Maurizot, J.C. 2005. New synchrotron radiation circular dichroism end-station on DISCO beamline at SOLEIL synchrotron for biomolecular analysis. *Biochim. Biophys. Acta* 1724: 425-431.

CHROMOSOMAL LOCATION

Genetic locus: GALM (human) mapping to 2p22.1.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PRODUCT

Galactose Mutarotase (h): 293 Lysate represents a lysate of human Galactose Mutarotase transfected 293 cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

Galactose Mutarotase (h): 293 Lysate is suitable as a Western Blotting positive control for human reactive Galactose Mutarotase antibodies. Recommended use: 10-20 μ l per lane.


Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

Galactose Mutarotase (H-7): sc-166304 is recommended as a positive control antibody for Western Blot analysis of enhanced human Galactose Mutarotase expression in Galactose Mutarotase transfected 293 cells (starting dilution 1:100, dilution range 1:100-1:1,000).


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



Galactose Mutarotase (H-7): sc-166304. Western blot analysis of Galactose Mutarotase expression in non-transfected: sc-110760 (A) and human Galactose Mutarotase transfected: sc-112214 (B) 293 whole cell lysates.



Galactose Mutarotase (G-7): sc-376233. Western blot analysis of Galactose Mutarotase expression in non-transfected: sc-110760 (A) and human Galactose Mutarotase transfected: sc-112214 (B) 293 whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.