

MOAP1 (h): 293 Lysate: sc-112238

BACKGROUND

MOAP1 (modulator of apoptosis 1) is a 352 amino acid protein encoded by the human gene MOAP1. MOAP1 belongs to the PNMA family and contains one BH3-like domain and one RASSF1-binding domain. It is required for death receptor-dependent apoptosis. When MOAP1 is associated with RASSF1, it promotes a Bax conformational change and translocation to mitochondrial membranes in response to TNF and TNFSF10 stimulation. MOAP1 is a homodimer and, under normal circumstances, is held in an inactive conformation by an intramolecular interaction. Binding to RASSF1 isoform A (RASSF1A) relieves this inhibitory interaction and allows further binding to Bax. MOAP1 will also bind to Bcl-2 and Bcl-x.

REFERENCES

1. Tan, K.O., Tan, K.M., Chan, S.L., Yee, K.S., Bevort, M., Ang, K.C. and Yu, V.C. 2001. MAP-1, a novel proapoptotic protein containing a BH3-like motif that associates with Bax through its Bcl-2 homology domains. *J. Biol. Chem.* 276: 2802-2807.
2. Tan, K.O., Fu, N.Y., Sukumaran, S.K., Chan, S.L., Kang, J.H., Poon, K.L., Chen, B.S. and Yu, V.C. 2005. MAP-1 is a mitochondrial effector of Bax. *Proc. Natl. Acad. Sci. USA* 102: 14623-14628.
3. Baksh, S., Tommasi, S., Fenton, S., Yu, V.C., Martins, L.M., Pfeifer, G.P., Latif, F., Downward, J. and Neel, B.G. 2005. The tumor suppressor RASSF1A and MAP-1 link death receptor signaling to Bax conformational change and cell death. *Mol. Cell* 18: 637-650.
4. Tretyakova, I., Zolotukhin, A.S., Tan, W., Bear, J., Propst, F., Ruthel, G. and Felber, B.K. 2005. Nuclear export factor family protein participates in cytoplasmic mRNA trafficking. *J. Biol. Chem.* 280: 31981-31990.
5. Schüller, M., Jenne, D. and Voltz, R. 2005. The human PNMA family: novel neuronal proteins implicated in paraneoplastic neurological disease. *J. Neuroimmunol.* 169: 172-176.
6. Vos, M.D., Dallol, A., Eckfeld, K., Allen, N.P., Donniger, H., Hesson, L.B., Calvisi, D., Latif, F. and Clark, G.J. 2006. The RASSF1A tumor suppressor activates Bax via MOAP1. *J. Biol. Chem.* 281: 4557-4563.
7. Fu, N.Y., Sukumaran, S.K. and Yu, V.C. 2007. Inhibition of ubiquitin-mediated degradation of MOAP1 by apoptotic stimuli promotes Bax function in mitochondria. *Proc. Natl. Acad. Sci. USA* 104: 10051-10056.

CHROMOSOMAL LOCATION

Genetic locus: MOAP1 (human) mapping to 14q32.12.

PRODUCT

MOAP1 (h): 293 Lysate represents a lysate of human MOAP1 transfected 293 cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

MOAP1 (h): 293 Lysate is suitable as a Western Blotting positive control for human reactive MOAP1 antibodies. Recommended use: 10-20 µl per lane.

Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.


MOAP1 (A-11): sc-271338 is recommended as a positive control antibody for Western Blot analysis of enhanced human MOAP1 expression in MOAP1 transfected 293 cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:

- 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



MOAP1 (A-11): sc-271338. Western blot analysis of MOAP1 expression in non-transfected: sc-110760 (A) and human MOAP1 transfected: sc-112238 (B) 293 whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.