

SOD-2 (h): 293 Lysate: sc-112776

BACKGROUND

The superoxide dismutase family is composed of three metalloenzymes (SOD-1, SOD-2 and SOD-3) that catalyze the oxido-reduction of reactive oxygen species (ROS), such as superoxide anion. The SOD-2 precursor is a 222 amino acid protein that is encoded by nuclear chromatin, synthesized in the cytosol and imported posttranslationally into the mitochondrial matrix. Unlike SOD-1, which is a homodimeric cytosolic Cu-Zn enzyme, SOD-2 is a homotetrameric manganese enzyme (also known as MnSOD) that functions in the mitochondrion. ROSs are implicated in a wide range of degenerative processes, including Alzheimer's disease, Parkinson's disease and ischemic heart disease. Homozygous mutant mice, which lack SOD-2, exhibit dilated cardiomyopathy, accumulation of lipid in liver and skeletal muscle, metabolic acidosis, oxidative DNA damage and respiratory chain deficiencies in heart and skeletal muscle. Polymorphisms in the SOD-2 gene have also been implicated in nonfamilial, idiopathic, dilated cardiomyopathy in humans.

REFERENCES

1. Wispé, J.R., Clark, J.C., Burhans, M.S., Kropp, K.E., Korfhagen, T.R. and Whitsett, J.A. 1989. Synthesis and processing of the precursor for human manganese-superoxide dismutase. *Biochem. Biophys. Acta* 994: 30-36.
2. Nishi, H., Koga, Y., Koyanagi, T., Harada, H., Imaizumi, T., Toshima, H., Sasazuki, T. and Kimura, A. 1995. DNA typing of HLA class II genes in Japanese patients with dilated cardiomyopathy. *J. Mol. Cell. Cardiol.* 27: 2385-2392.
3. Li, Y., Huang, T.T., Carlson, E.J., Melov, S., Ursell, P.C., Olson, J.L., Noble, L.J., Yoshimura, M.P., Berger, C., Chan, P.H., Wallace, D.C. and Epstein, C.J. 1995. Dilated cardiomyopathy and neonatal lethality in mutant mice lacking manganese superoxide dismutase. *Nat. Genet.* 11: 376-381.
4. Borgstahl, G.E., Parge, H.E., Hickey, M.J., Johnson, M.J., Boissinot, M., Hallewell, R.A., Lepock, J.R., Cabelli, D.E. and Tainer, J.A. 1996. Human mitochondrial manganese superoxide dismutase polymorphic variant Ile58Thr reduces activity by destabilizing the tetrameric interface. *Biochemistry* 35: 4287-4297.
5. Hsieh, Y., Guan, Y., Tu, C., Bratt, P.J., Angerhofer, A., Lepock, J.R., Hickey, M.J., Tainer, J.A., Nick, H.S. and Silverman, D.N. 1998. Probing the active site of human manganese superoxide dismutase: the role of glutamine 143. *Biochemistry* 37: 4731-4739.
6. Melov, S., Schneider, J.A., Day, B.J., Hinerfeld, D., Coskun, P., Mirra, S.S., Crapo, J.D. and Wallace, D.C. 1998. A novel neurological phenotype in mice lacking mitochondrial manganese superoxide dismutase. *Nat. Genet.* 18: 159-163.
7. Melov, S., Coskun, P., Patel, M., Tuinstra, R., Cottrell, B., Jun, A.S., Zastawny, T.H., Dizdaroglu, M., Goodman, S.I., Huang, T.T., Miziorko, H., Epstein, C.J. and Wallace, D.C. 1999. Mitochondrial disease in superoxide dismutase 2 mutant mice. *Proc. Natl. Acad. Sci. USA* 96: 846-851.
8. Hiroi, S., Harada, H., Nishi, H., Satoh, M., Nagai, R. and Kimura, A. 1999. Polymorphisms in the SOD-2 and HLA-DRB1 genes are associated with non-familial idiopathic dilated cardiomyopathy in Japanese. *Biochem. Biophys. Res. Commun.* 261: 332-339.

CHROMOSOMAL LOCATION

Genetic locus: SOD2 (human) mapping to 6q25.3.

PRODUCT

SOD-2 (h): 293 Lysate represents a lysate of human SOD-2 transfected 293 cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

SOD-2 (h): 293 Lysate is suitable as a Western Blotting positive control for human reactive SOD-2 antibodies. Recommended use: 10-20 µl per lane.

Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.