

SIRT1 (h): 293T Lysate: sc-113797

BACKGROUND

The silent information regulator (SIR) family of genes are highly-conserved from prokaryotes to eukaryotes and are involved in diverse processes, including transcriptional regulation, cell cycle progression, DNA-damage repair and aging. In *S. cerevisiae*, Sir2p deacetylates histones in a NAD-dependent manner, which regulates silencing at the telomeric, rDNA and silent mating-type loci. Sir2p is the founding member of a large family, designated sirtuins, which contain a conserved catalytic domain. The human homologs, which include SIRT1-7, are divided into four main branches: SIRT1-3 are class I, SIRT4 is class II, SIRT5 is class III and SIRT6-7 are class IV. SIRT1 has the closest homology to the yeast Sir2p and is widely expressed in fetal and adult tissues. SIRT1 is highly expressed in heart, brain and skeletal muscle, with low expression in lung and placenta. SIRT1 regulates the p53-dependent DNA damage response pathway by binding to and deacetylating p53, specifically at Lys 382.

REFERENCES

1. Frye, R.A. 1999. Characterization of five human cDNAs with homology to the yeast SIR2 gene: Sir2-like proteins (sirtuins) metabolize NAD and may have protein ADP-ribosyltransferase activity. *Biochem. Biophys. Res. Commun.* 260: 273-279.
2. Afshar, G. and Murnane, J.P. 1999. Characterization of a human gene with sequence homology to *Saccharomyces cerevisiae* Sir2. *Gene* 234: 161-168.
3. Sherman, J.M., Stone, E.M., Freeman-Cook, L.L., Brachmann, C.B., Boeke, J.D. and Pillus, L. 1999. The conserved core of a human SIR2 homolog functions in yeast silencing. *Mol. Biol. Cell* 10: 3045-3059.
4. Frye, R.A. 2000. Phylogenetic classification of prokaryotic and eukaryotic Sir2-like proteins. *Biochem. Biophys. Res. Commun.* 273: 793-798.
5. Grozinger, C.M., Chao, E.D., Blackwell, H.E., Moazed, D. and Schreiber, S.L. 2001. Identification of a class of small molecule inhibitors of the sirtuin family of NAD-dependent deacetylases by phenotypic screening. *J. Biol. Chem.* 276: 38837-38843.
6. Defossez, P.A., Lin, S.J. and McNabb, D.S. 2001. Sound silencing: the Sir2 protein and cellular senescence. *Bioessays* 23: 327-332.
7. Vaziri, H., Dessain, S.K., Ng Eaton, E., Imai, S.I., Frye, R.A., Pandita, T.K., Guarente, L. and Weinberg, R.A. 2001. hSIR2 (SIRT1) functions as a NAD-dependent p53 deacetylase. *Cell* 107: 149-159.

CHROMOSOMAL LOCATION

Genetic locus: SIRT1 (human) mapping to 10q21.3.

PRODUCT

SIRT1 (h): 293T Lysate represents a lysate of human SIRT1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

SIRT1 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive SIRT1 antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

SIRT1 (B-7): sc-74465 is recommended as a positive control antibody for Western Blot analysis of enhanced human SIRT1 expression in SIRT1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



SIRT1 (B-7): sc-74465. Western blot analysis of SIRT1 expression in non-transfected: sc-117752 (A) and human SIRT1 transfected: sc-113797 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.