

CD32 (h3): 293T Lysate: sc-113838

BACKGROUND

CD32 (also designated FcγRII) is a low affinity receptor for the Fc fragment of aggregated IgG. CD32 is responsible for the clearance of immunocomplexes by macrophages and also plays an important role in the regulation of antibody production by B cells. IgG can noncooperatively bind either one or two highly-glycosylated CD32 molecules, and this binding delivers a negative signal for B cells. CD32 exists as several isoforms that are produced by alternative splicing of three distinct genes; A, B and C. These isoforms are designated FcγRIIA, FcγRIIB1, FcγRIIB3 and FcγRIIC. All isoforms are present on monocytes, placental trophoblasts and endothelial cells. In addition, the FcγRIIB forms are present on B lymphocytes, and the FcγRIIA and FcγRIIC forms are found on neutrophils.

REFERENCES

1. Bijstervosch, M.K. and Klaus, G.G. 1985. Crosslinking of surface immunoglobulin and Fc receptors on B lymphocytes inhibits stimulation of inositol phospholipid breakdown via the antigen receptors. *J. Exp. Med.* 162: 1825-1836.
2. Huizinga, T.W.J., Kerst, M., Nuyens, J.H., Vlug, A., von dem Borne, A.E., Roos, D. and Tetteroo, P.A. 1989. Binding characteristics of dimeric IgG subclass complexes to human neutrophils. *J. Immunol.* 142: 2365-2369.
3. Stuart, S.G., Simister, N.E., Clarkson, S.B., Kacinski, B.M., Shapiro, M. and Mellman, I. 1989. Human IgG Fc receptor (hFcRII; CD32) exists as multiple isoforms in macrophages, lymphocytes and IgG-transporting placental epithelium. *EMBO J.* 8: 3657-3666.
4. van de Winkel, J.G.J. and Anderson, C.L. 1991. CD32 cluster workshop report. In Schlossman SF, Boumsell, L., Gilks, W., Harlan, J.M., Kishimoto, T., Morimoto, C., Ritz, J., Shaw, S., Silverstein, R.L., Springer, T.A., Tedder, T.F. and Todd, R.F.(eds) *Leucocyte Typing V*, Vol 1, Oxford University Press, Oxford, New York and Tokyo, p 823.
5. Raveth, J.V. and Kinet, J.P. 1991. Fc receptors. *Annu. Rev. Immunol.* 9: 457-492.
6. Barclay, A.N., Beyers, A.D., Birkeland, M.L., Brown, S.J., Somoza, C. and Williams, A.F. 1993. *The Leukocyte Antigen Facts Book*. London: Academic Press, 170-172.
7. Budde, P., Weinrich, V., Sondermann, P., Bewarder, N., Kilian, A., Schulzeck, O. and Frey, J. 1995. Specificity of CD32 mAb for FcγRIIA, FcγRIIB1, and FcγRIIB2 expressed in transfected mouse B cells and BHK-21 cells: section report. In Schlossman S.F, Boumsell, L., Gilks, W., Harlan, J.M., Kishimoto, T., Morimoto, C., Ritz, J., Shaw, S., Silverstein, R.L., Springer, T.A., Tedder, T.F. and Todd, R.F.(eds) *Leucocyte Typing V*, Vol 1, Oxford University Press, Oxford, New York and Tokyo, p 828.
8. Sondermann, P., Jacob, U., Kutscher, C. and Frey, J. 1999. Characterization and crystallization of soluble human FcγRII (CD32) isoforms produced in insect cells. *Biochemistry* 38: 8469-8477.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: FCGR2A/FCGR2B/FCGR2C (human) mapping to 1q23.3.

PRODUCT

CD32 (h3): 293T Lysate represents a lysate of human CD32 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

CD32 (h3): 293T Lysate is suitable as a Western Blotting positive control for human reactive CD32 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


CD32-A/B/C (B-4) HRP: sc-166711 HRP is recommended as a positive control antibody for Western Blot analysis of enhanced human CD32 expression in CD32 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS


To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



CD32-A/B/C (B-4) HRP: sc-166711 HRP. Direct western blot analysis of CD32 expression in non-transfected: sc-117752 (A) and human CD32 transfected: sc-113838 (B) 293T whole cell lysates.



CD32-B (F-4): sc-365864. Western blot analysis of CD32 expression in non-transfected: sc-117752 (A) and human CD32 transfected: sc-113838 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.