

connexin 30.3 (h): 293T Lysate: sc-113859

BACKGROUND

The connexin family of proteins form hexameric complexes called "connexons" that facilitate movement of low molecular weight proteins between cells via gap junctions. Connexin proteins share a common topology of four transmembrane α -helical domains, two extracellular loops, a cytoplasmic loop and cytoplasmic N- and C-termini. Many of the key functional differences arise from specific amino-acid substitutions in the most highly conserved domains: the transmembrane and extracellular regions. Connexin 30.3, also known as GJB4 (gap junction β -4 protein), CX30.3 or EKV, is expressed in the epidermis (upper spinous and granular layers) and is believed to play a role in keratinocyte intercellular communication. Mutations in the gene encoding connexin 30.3 can result in erythrokeratoderma variabilis (EVK), a condition characterized by localized or generalized hyperkeratosis and random, transient erythematous patches.

REFERENCES

- Manjunath, C.K., et al. 1987. Human cardiac gap junctions: isolation, ultrastructure, and protein composition. *J. Mol. Cell. Cardiol.* 19: 131-134.
- Macari, F., et al. 2000. Mutation in the gene for connexin 30.3 in a family with erythrokeratoderma variabilis. *Am. J. Hum. Genet.* 67: 1296-1301.
- Harris, A.L. 2001. Emerging issues of connexin channels: biophysics fills the gap. *Q. Rev. Biophys.* 34: 325-472.
- Plantard, L., et al. 2003. Molecular interaction of connexin 30.3 and connexin 31 suggests a dominant-negative mechanism associated with erythrokeratoderma variabilis. *Hum. Mol. Genet.* 12: 3287-3294.
- Richard, G., et al. 2003. Genetic heterogeneity in erythrokeratoderma variabilis: novel mutations in the connexin gene GJB4 (CX30.3) and genotype-phenotype correlations. *J. Invest. Dermatol.* 120: 601-609.
- Arita, K., et al. 2003. Erythro-keratoderma variabilis without connexin 31 or connexin 30.3 gene mutation: immunohistological, ultrastructural and genetic studies. *Acta Derm. Venereol.* 83: 266-270.
- Fonseca, P.C., et al. 2004. Characterization of connexin 30.3 and 43 in thymocytes. *Immunol. Lett.* 94: 65-75.
- Tamaki, Y., et al. 2006. A case of erythrokeratoderma variabilis: loosened gap junctions in the acanthotic epidermis. *J. Dermatol.* 33: 419-423.
- Zheng-Fischhöfer, Q., et al. 2007. Characterization of connexin 30.3-deficient mice suggests a possible role of connexin 30.3 in olfaction. *Eur. J. Cell Biol.* 86: 683-700.

CHROMOSOMAL LOCATION

Genetic locus: GJB4 (human) mapping to 1p34.3.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

connexin 30.3 (h): 293T Lysate represents a lysate of human connexin 30.3 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

connexin 30.3 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive connexin 30.3 antibodies. Recommended use: 10-20 μ l per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

connexin 30.3 (F6): sc-81802 is recommended as a positive control antibody for Western Blot analysis of enhanced human connexin 30.3 expression in connexin 30.3 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



connexin 30.3 (F6): sc-81802. Western blot analysis of connexin 30.3 expression in non-transfected: sc-117752 (A) and human connexin 30.3 transfected: sc-113859 (B) 293T whole cell lysates.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.