



Importin-11 (h): 293T Lysate: sc-114980

BACKGROUND

The Importin complex consists of Importin- α and Importin- β proteins which assist in the transport of arginine- or serine-rich (SR) peptides across the nucleus. Importin-11, also known as IPO11 or RanBP11, is a 975 amino acid protein that contains one importin N-terminal domain and 15 HEAT repeats and belongs to the Importin- β family. Localized to both the nucleus and the cytoplasm, Importin-11 interacts with UBE2E3 and plays a role in nuclear protein import, specifically functioning as a nuclear transport receptor that mediates the docking of the importin complex to the nuclear pore complex (NPC). The gene encoding Importin-11 maps to human chromosome 5, which contains 181 million base pairs and comprises nearly 6% of the human genome. Deletion of the p arm of chromosome 5q12.1 leads to Cri du chat syndrome, while deletion of the q arm of chromosome 5q12.1 altogether is common in therapy-related acute myelogenous leukemias and myelodysplastic syndrome.

REFERENCES

1. Plafker, S.M. and Macara, I.G. 2000. Importin-11, a nuclear import receptor for the ubiquitin-conjugating enzyme, UbcM2. *EMBO J.* 19: 5502-5513.
2. Plafker, S.M. and Macara, I.G. 2002. Ribosomal protein L12 uses a distinct nuclear import pathway mediated by Importin-11. *Mol. Cell. Biol.* 22: 1266-1275.
3. Wang, W., Yang, Y., Li, L. and Shi, Y. 2003. Synleucin, a novel leucine-rich repeat protein that increases the intensity of pleiotropic cytokine responses. *Biochem. Biophys. Res. Commun.* 305: 981-988.
4. Plafker, S.M., Plafker, K.S., Weissman, A.M. and Macara, I.G. 2004. Ubiquitin charging of human class III ubiquitin-conjugating enzymes triggers their nuclear import. *J. Cell Biol.* 167: 649-659.
5. Zhang, X.D. and Matunis, M.J. 2005. Ub in charge: regulating E2 enzyme nuclear import. *Nat. Cell Biol.* 7: 12-14.
6. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 610889. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: IPO11 (human) mapping to 5q12.1.

PRODUCT

Importin-11 (h): 293T Lysate represents a lysate of human Importin-11 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

APPLICATIONS

Importin-11 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive Importin-11 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.