

Mx1 (h): 293T Lysate: sc-115203

BACKGROUND

The Dynamin family of microtubule-associated proteins function as GTPases that are involved in microtubule bundling and endocytosis. In mice, Mx2 (myxovirus resistance protein 2) and Mx1 (myxovirus resistance protein 1) are members of the Dynamin family that are involved in the immune response to viral infections. Localized to the cytoplasm, Mx2 contains one GED domain and is expressed in response to viral infection or treatment by IFN- α /IFN- β . Once expression is induced, Mx2 accumulates in the cytoplasm and inhibits the replication of vesicular stomatitis virus (VSV), thereby conferring resistance to VSV infection. Unlike Mx2, Mx1 is localized to the nucleus where, upon induction by IFN- α /IFN- β , it provides selective resistance to infection by the highly lethal H5N1 influenza virus. In humans, MxA and MxB function in a similar manner to Mx1 and Mx2, conferring resistance to specific target viruses.

REFERENCES

1. Lindenmann, J. 1964. Inheritance of resistance to influenza virus in mice. *Proc. Soc. Exp. Biol. Med.* 116: 506-509.
2. Staeheli, P., Haller, O., Boll, W., Lindenmann, J. and Weissmann, C. 1986. Mx protein: constitutive expression in 3T3 cells transformed with cloned Mx cDNA confers selective resistance to influenza virus. *Cell* 44: 147-158.
3. Hug, H., Costas, M., Staeheli, P., Aebi, M. and Weissmann, C. 1988. Organization of the murine Mx gene and characterization of its interferon- and virus-inducible promoter. *Mol. Cell. Biol.* 8: 3065-3079.
4. Staeheli, P. and Sutcliffe, J.G. 1988. Identification of a second interferon-regulated murine Mx gene. *Mol. Cell. Biol.* 8: 4524-4528.
5. Zürcher, T., Pavlovic, J. and Staeheli, P. 1992. Mouse Mx2 protein inhibits vesicular stomatitis virus but not influenza virus. *Virology* 187: 796-800.
6. Jin, H.K., Takada, A., Kon, Y., Haller, O. and Watanabe, T. 1999. Identification of the murine Mx2 gene: interferon-induced expression of the Mx2 protein from the feral mouse gene confers resistance to vesicular stomatitis virus. *J. Virol.* 73: 4925-4930.
7. Lee, S.H. and Vidal, S.M. 2002. Functional diversity of Mx proteins: variations on a theme of host resistance to infection. *Genome Res.* 12: 527-530.
8. Salomon, R., Staeheli, P., Kochs, G., Yen, H.L., Franks, J., Reh, J.E., Webster, R.G. and Hoffmann, E. 2007. Mx1 gene protects mice against the highly lethal human H5N1 influenza virus. *Cell Cycle* 6: 2417-2421.
9. Tumpey, T.M., Szretter, K.J., Van Hoeven, N., Katz, J.M., Kochs, G., Haller, O., García-Sastre, A. and Staeheli, P. 2007. The Mx1 gene protects mice against the pandemic 1918 and highly lethal human H5N1 influenza viruses. *J. Virol.* 81: 10818-10821.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: MX1 (human) mapping to 21q22.3.

PRODUCT

Mx1 (h): 293T Lysate represents a lysate of human Mx1 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

Mx1 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive Mx1 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


Mx1/2/3 (C-1): sc-166412 is recommended as a positive control antibody for Western Blot analysis of enhanced human Mx1 expression in Mx1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS


To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



Mx1/2/3 (C-1): sc-166412. Western blot analysis of Mx1 expression in non-transfected: sc-117752 (A) and human Mx1 transfected: sc-115203 (B) 293T whole cell lysates.



Mx1 (E-5): sc-271024. Near-infrared western blot analysis of Mx1 expression in non-transfected: sc-117752 (A) and human Mx1 transfected: sc-115203 (B) 293T whole cell lysates. Blocked with UltraCruz® Blocking Reagent: sc-516214. Detection reagent used: m-IgG κ BP-CFL 680: sc-516180.

RESEARCH USE

For research use only, not for use in diagnostic procedures.