

TULP3 (h): 293T Lysate: sc-115216

BACKGROUND

Mutations in the mouse *Tub* gene gradually lead to obesity, strongly resembling the late-onset obesity observed in the human population. In addition to excessive deposition of adipose tissue, mice with the *Tub* phenotype also suffer retinal degeneration and neurosensory hearing loss. A human homolog of the *Tub* gene has been identified, as have three related proteins, called tubby-like protein 1 (TULP1), TULP2 and TULP3. When compared to TULP1 and TULP2, TULP3 has a wider tissue expression and is phylogenetically more similar to *Tub* than either TULP1 or TULP2. TULP1, expressed specifically in the retina, maps to the chromosomal region known to be involved in retinitis pigmentosa, while TULP2 maps within the minimal interval for the rod-cone dystrophy. TULP3 maps to human chromosome 12p13.33, and shares 69% homology to mouse TULP3. Human RNA from testis, ovary, thyroid and spinal cord contain highly detectable levels of TULP3 transcripts. In the retina, TULP3 is expressed specifically in the inner nuclear layer and ganglion cell layer. TULP1, TULP2 and TULP3 may comprise a unique family of bipartite transcription factors.

REFERENCES

1. Ohlemiller, K.K., et al. 1995. Cochlear and retinal degeneration in the tubby mouse. *Neuroreport* 6: 845-849.
2. Noben-Trauth, K., et al. 1996. A candidate gene for the mouse mutation tubby. *Nature* 380: 534-538.
3. North, M.A., et al. 1997. Molecular characterization of TUB, TULP1, and TULP2, members of the novel tubby gene family and their possible relation to ocular diseases. *Proc. Natl. Acad. Sci. USA* 94: 3128-3133.
4. Gu, S., et al. 1998. Tubby-like protein-1 mutations in autosomal recessive retinitis pigmentosa. *Lancet* 351: 1103-1104.
5. Hagstrom, S.A., et al. 1998. Recessive mutations in the gene encoding the tubby-like protein TULP1 in patients with retinitis pigmentosa. *Nat. Genet.* 18: 174-176.
6. Banerjee, P., et al. 1998. TULP1 mutation in two extended Dominican kindreds with autosomal recessive retinitis pigmentosa. *Nat. Genet.* 18: 177-179.
7. Nishina, P.M., et al. 1998. Molecular characterization of a novel tubby gene family member, TULP3, in mouse and humans. *Genomics* 54: 215-220.
8. Ikeda, S., et al. 1999. Cells-specific expression of tubby gene family members (Tub, Tulp1, 2, and 3) in the retina. *Invest. Ophthalmol. Vis. Sci.* 40: 2706-2712.
9. Boggon, T.J., et al. 1999. Implication of tubby proteins as transcription factors by structure-based functional analysis. *Science* 286: 2119-2125.

CHROMOSOMAL LOCATION

Genetic locus: TULP3 (human) mapping to 12p13.33.

PRODUCT

TULP3 (h): 293T Lysate represents a lysate of human TULP3 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

TULP3 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive TULP3 antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

TULP3 (G-4): sc-514733 is recommended as a positive control antibody for Western Blot analysis of enhanced human TULP3 expression in TULP3 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



TULP3 (G-4): sc-514733. Western blot analysis of TULP3 expression in non-transfected: sc-117752 (A) and human TULP3 transfected: sc-115216 (B) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.