

dUTPase (h2): 293T Lysate: sc-115300

BACKGROUND

dUTPase (deoxyuridine 5'-triphosphate nucleotidohydrolase), also known as DUT or dUTP pyrophosphatase, is a preventive DNA repair enzyme that functions in nucleotide metabolism. dUTPase is expressed in a variety of tissues and, depending on the isoform (DUT-N or DUT-M), localizes to the nucleus or the mitochondrion. The nuclear isoform, DUT-N, is the most abundant of the two isoforms. dUTPase, in the presence of magnesium ions, is responsible for hydrolyzing dUTP to dUMP and diphosphate. This reaction is important for keeping the intracellular dUTP concentration low so that uracil does not become incorporated into DNA. Extensive incorporation of uracil into DNA can ultimately lead to cell death. This suggests that dUTPase is essential for cell viability, further implying that dUTPase is a potential target for anticancer therapy. In addition, dUMP, the product of the hydrolysis reaction, is a precursor of thymidine nucleotides which are essential for DNA replication.

REFERENCES

- Canman, C.E., Tang, H.Y., Normolle, D.P., Lawrence, T.S. and Maybaum, J. 1992. Variations in patterns of DNA damage induced in human colorectal tumor cells by 5-fluorodeoxyuridine: implications for mechanisms of resistance and cytotoxicity. *Proc. Natl. Acad. Sci. USA* 89: 10474-10478.
- Ladner, R.D., McNulty, D.E., Carr, S.A., Roberts, G.D. and Caradonna, S.J. 1996. Characterization of distinct nuclear and mitochondrial forms of human deoxyuridine triphosphate nucleotidohydrolase. *J. Biol. Chem.* 271: 7745-7751.
- Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 601266. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Jiang, Y.L., Chung, S., Krosky, D.J. and Stivers, J.T. 2006. Synthesis and high-throughput evaluation of triskelion uracil libraries for inhibition of human dUTPase and UNG2. *Bioorg. Med. Chem.* 14: 5666-5672.
- Samal, A., Schormann, N., Cook, W.J., DeLucas, L.J. and Chattopadhyay, D. 2007. Structures of vaccinia virus dUTPase and its nucleotide complexes. *Acta Crystallogr. D Biol. Crystallogr.* 63: 571-580.
- Varga, B., Barabás, O., Kovári, J., Tóth, J., Hunyadi-Gulyás, E., Klement, E., Medzihradsky, K.F., Tölgyesi, F., Fidy, J. and Vértessy, B.G. 2007. Active site closure facilitates juxtaposition of reactant atoms for initiation of catalysis by human dUTPase. *FEBS Lett.* 581: 4783-4788.
- Tóth, J., Varga, B., Kovács, M., Málnási-Csizmadia, A. and Vértessy, B.G. 2007. Kinetic mechanism of human dUTPase, an essential nucleotide pyrophosphatase enzyme. *J. Biol. Chem.* 282: 33572-33582.
- Kovári, J., Barabás, O., Varga, B., Békési, A., Tölgyesi, F., Fidy, J., Nagy, J. and Vértessy, B.G. 2007. Methylene substitution at the α - β bridging position within the phosphate chain of dUDP profoundly perturbs ligand accommodation into the dUTPase active site. *Proteins* 71: 308-319.
- Thymark, M., Johansson, E., Larsen, S. and Willemoës, M. 2008. Mutational analysis of the nucleotide binding site of *Escherichia coli* dCTP deaminase. *Arch. Biochem. Biophys.* 470: 20-26.

CHROMOSOMAL LOCATION

Genetic locus: DUT (human) mapping to 15q21.1.

PRODUCT

dUTPase (h2): 293T Lysate represents a lysate of human dUTPase transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

dUTPase (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive dUTPase antibodies. Recommended use: 10-20 μ l per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

dUTPase (G-6): sc-166838 is recommended as a positive control antibody for Western Blot analysis of enhanced human dUTPase expression in dUTPase transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



dUTPase (G-6): sc-166838. Western blot analysis of dUTPase expression in non-transfected: sc-117752 (A) and human dUTPase transfected: sc-115300 (B) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.