



MSY2 (h2): 293T Lysate: sc-116103

BACKGROUND

MSY2 and YB-2 (MSY3,4) belong to the Y-box family of multifunctional proteins that regulate both transcription and translation. Y-box proteins interact with a wide variety of nucleic acid structures to act as transcription factors and mRNA masking proteins. The modular structure of Y-box proteins includes a highly conserved N-terminal cold-shock domain (CSD, equivalent to the bacterial cold-shock proteins) and four basic C-terminal domains containing arginine clusters and aromatic residues. MSY2 is expressed in testis and ovary where it may repress translation of parental mRNA. The gene encoding human MSY2 maps to chromosome 17p13.1. YB-2 (MSY3,4 in mouse) is also known as DNA binding protein A and is highly expressed in the testis, heart and muscle. MSY2 and YB-2 bind to the consensus sequence 5'-UCCAUCA-3' contained in the Y-box element.

REFERENCES

1. Tekur, S., et al. 1999. Contrin, the human homologue of a germ-cell Y-box-binding protein: cloning, expression, and chromosomal localization. *J. Androl.* 20: 135-144.
2. Okamoto, T., et al. 2000. Direct interaction of p53 with the Y-box binding protein, YB-1: a mechanism for regulation of human gene expression. *Oncogene* 54: 6194-6202.
3. Davies, H.G., et al. 2000. A sequence-specific RNA binding complex expressed in murine germ cells contains MSY2 and MSY4. *Dev. Biol.* 221: 87-100.
4. Levenson, V.V., et al. 2000. Pleiotropic resistance to DNA-interactive drugs is associated with increased expression of genes involved in DNA replication, repair, and stress response. *Cancer Res.* 18: 5027-5030.
5. Mastrangelo, M.A. and Kleene, K.C. 2000. Developmental expression of Y-box protein 1 mRNA and alternatively spliced Y-box protein 3 mRNAs in spermatogenic cells in mice. *Mol. Hum. Reprod.* 6: 779-788.
6. Diamond, P., et al. 2001. Cold shock domain factors activate the granulocyte-macrophage colony-stimulating factor promoter in stimulated Jurkat T cells. *J. Biol. Chem.* 11: 7943-7951.
7. Chansky, H.A., et al. 2001. Oncogenic TLS/ERG and EWS/Fli-1 fusion proteins inhibit RNA splicing mediated by YB-1 protein. *Cancer Res.* 9: 3586-3590.

CHROMOSOMAL LOCATION

Genetic locus: YBX2 (human) mapping to 17p13.1.

PRODUCT

MSY2 (h2): 293T Lysate represents a lysate of human MSY2 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

MSY2 (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive MSY2 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.