

SREBP-1 (h): 293T Lysate: sc-116164

BACKGROUND

The low density lipoprotein (LDL) receptor mediates the endocytic uptake of cholesterol-carrying lipoproteins, thereby controlling cholesterol levels in cells and plasma. Transcription of the LDL receptor gene is controlled by a 10 base pair sequence in the 5' flanking region, designated sterol regulatory element 1 (SRE-1). When cellular sterol stores are depleted, the element is activated, the gene is transcribed and the cellular uptake of LDL increases. A set of SRE-binding proteins (SREBPs) have been identified, including two basic helix-loop-helix leucine zipper (bHLH-Zip) transcription factors, designated SREBP-1 and SREBP-2. SREBP-1 (also designated ADD1, for adipocyte determination and differentiation factor) is synthesized as a precursor that is attached to the nuclear envelope and endoplasmic reticulum. In sterol-depleted cells, the membrane-bound precursor is cleaved to generate a soluble NH₂-terminal fragment that translocates to the nucleus to activate transcription. Sterols inhibit the cleavage of SREBP-1.

REFERENCES

1. Brown, M.S. and Goldstein, J.L. 1986. A receptor-mediated pathway for cholesterol homeostasis. *Science* 232: 34-47.
2. Smith, J.R., Osborne, T.F., Goldstein, J.L. and Brown, M.S. 1990. Identification of nucleotides responsible for enhancer activity of sterol regulatory element in low density lipoprotein receptor gene. *J. Biol. Chem.* 265: 2306-2310.
3. Goldstein, J.L. and Brown, M.S. 1990. Regulation of the mevalonate pathway. *Nature* 343: 425-430.
4. Briggs, M.R., Yokoyama, C., Wang, X., Brown, M.S. and Goldstein, J.L. 1993. Nuclear protein that binds sterol regulatory element of LDL receptor promoter. I. Identification of the protein and delineation of its target nucleotide sequence. *J. Biol. Chem.* 268: 14490-14496.
5. Wang, X., Briggs, M.R., Hua, X., Yokoyama, C., Goldstein, J.L. and Brown, M.S. 1993. Nuclear protein that binds sterol regulatory element of LDL receptor promoter. II. Purification and characterization. *J. Biol. Chem.* 268: 14497-14504.
6. Hua, X., Yokoyama, C., Wu, J., Briggs, M.R., Brown, M.S., Goldstein, J.L. and Wang, X. 1993. SREBP-2, a second basic-helix-loop-helix-leucine zipper protein that stimulates transcription by binding to a sterol regulatory element. *Proc. Natl. Acad. Sci. USA* 90: 11603-11607.
7. Wang, X., Sato, R., Brown, M.S., Hua, X. and Goldstein, J.L. 1994. SREBP-1, a membrane-bound transcription factor released by sterol-regulated proteolysis. *Cell* 77: 53-62.
8. Rebouissou, S., Imbeaud, S., Balabaud, C., Boulanger, V., Bertrand-Michel, J., Tercé, F., Auffray, C., Bioulac-Sage, P. and Zucman-Rossi, J. 2007. HNF-1 α inactivation promotes lipogenesis in human hepatocellular adenoma independently of SREBP-1 and ChREBP activation. *J. Biol. Chem.* 282: 14437-14446.
9. Harrison, W.J., Bull, J.J., Seltmann, H., Zouboulis, C.C. and Philpott, M.P. 2007. Expression of lipogenic factors galectin-12, resistin, SREBP-1, and SCD in human sebaceous glands and cultured sebocytes. *J. Invest. Dermatol.* 127: 1309-1317.

CHROMOSOMAL LOCATION

Genetic locus: SREBF1 (human) mapping to 17p11.2.

PRODUCT

SREBP-1 (h): 293T Lysate represents a lysate of human SREBP-1 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

SREBP-1 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive SREBP-1 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


SREBP-1 (E-4): sc-17755 is recommended as a positive control antibody for Western Blot analysis of enhanced human SREBP-1 expression in SREBP-1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS


To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



SREBP-1 (E-4): sc-17755. Western blot analysis of SREBP-1 expression in non-transfected: sc-117752 (A) and human SREBP-1 transfected: sc-116164 (B) 293T whole cell lysates.



SREBP-1 (SPM449): sc-56483. Western blot analysis of SREBP-1 expression in non-transfected: sc-117752 (A) and human SREBP-1 transfected: sc-116164 (B) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.