

TIF1 β (h): 293T Lysate: sc-116358

BACKGROUND

TIF1 β , for transcriptional intermediary factor 1 β , also designated KAP1 (for KRAB-associated protein 1), TF1 β and TRIM28 (for tripartite motif-containing 28), is a member of the tripartite motif family characterized by three zinc-binding domains, a RING finger, B-boxes and a coiled-coil domain. Like TIF1 α , TIF1 β contains both a Cys/His PHD (plant homeodomain) finger and bromodomain that form a cooperative unit required for transcriptional repression. TIF1 β mediates transcriptional control by interaction with the Kruppel-associated box (KRAB) repression domain found in many transcription factors and by binding DNA through its zinc finger. The human TIF1 β gene maps to human chromosome 19q13.4 and encodes an 835 amino acid nuclear protein.

REFERENCES

1. Friedman, J., Fredericks, W., Jensen, D., Speicher, D., Huang, X., Neilson, E. and Rauscher, F. 1996. KAP-1, a novel corepressor for the highly conserved KRAB repression domain. *Genes Dev.* 10: 2067-2078.
2. Moosmann, P., Georgiev, O., Le Douarin, B., Bourquin, J. and Schaffner, W. 1996. Transcriptional repression by RING finger protein TIF1 β that interacts with the KRAB repressor domain of KOX1. *Nucleic Acids Res.* 24: 4859-4867.
3. Schultz, D., Friedman, J. and Rauscher, F. 2001. Targeting histone deacetylase complexes via KRAB-zinc finger proteins: the PHD and bromodomains of KAP-1 form a cooperative unit that recruits a novel isoform of the Mi-2 α subunit of NuRD. *Genes Dev.* 15: 428-443.
4. Gebelein, B. and Urrutia, R. 2001. Sequence-specific transcriptional repression by KS1, a multiple-zinc-finger-Kruppel-associated box protein. *Mol. Cell. Biol.* 21: 928-939.
5. Online Mendelian Inheritance in Man, OMIM™. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 601742. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. LocusLink Report (LocusID: 10155). <http://www.ncbi.nlm.nih.gov/LocusLink/>

CHROMOSOMAL LOCATION

Genetic locus: TRIM28 (human) mapping to 19q13.4.

PRODUCT

TIF1 β (h): 293T Lysate represents a lysate of human TIF1 β transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

TIF1 β (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive TIF1 β antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.