

PNPase (h): 293T Lysate: sc-116464

BACKGROUND

Mitochondrial polynucleotide nucleotidyltransferase, also designated 3'-5' RNA exonuclease, OLD35, PNPase or PNPT1, is an evolutionally conserved protein in which the mouse protein shares 90% identity with the human version. PNPase participates in mRNA degradation and hydrolyzes single-stranded ribonucleotides in the 3' to 5' direction. PNPase forms homotrimers and is upregulated in response to interferon- β induction. The N-terminus of PNPase contains a putative mitochondrial targeting sequence; mutation analysis confirms that N-terminal sequences of PNPase target the protein to the mitochondria. Endogenous PNPase also co-localizes with a mitochondrial marker protein in HeLa cells.

REFERENCES

1. Bermúdez-Cruz, R.M., Fernández-Ramírez, F., Kameyama-Kawabe, L. and Montañez, C. 2005. Conserved domains in polynucleotide phosphorylase among eubacteria. *Biochimie* 87: 737-745.
2. Bollenbach, T.J., Lange, H., Gutierrez, R., Erhardt, M., Stern, D.B. and Gagliardi, D. 2005. RNR1, a 3'-5' exoribonuclease belonging to the RNR superfamily, catalyzes 3' maturation of chloroplast ribosomal RNAs in *Arabidopsis thaliana*. *Nucleic Acids Res.* 33: 2751-2563.
3. Oussenko, I.A., Abe, T., Ujiie, H., Muto, A. and Bechhofer, D.H. 2005. Participation of 3'-to-5' exoribonucleases in the turnover of *Bacillus subtilis* mRNA. *J. Bacteriol.* 187: 2758-2767.
4. Sarkar, D., Park, E.S., Emdad, L., Randolph, A., Valerie, K. and Fisher, P.B. 2005. Defining the domains of human polynucleotide phosphorylase (hPNPaseOLD-35) mediating cellular senescence. *Mol. Cell. Biol.* 25: 7333-7343.
5. Gewartowski, K., Tomecki, R., Muchowski, L., Dmochow Ska, A., Dzwonek, A., Malecki, M., Skurzak, H., Ostrowski, J. and Stepień, P.P. 2006. Up-regulation of human PNPase mRNA by β -interferon has no effect on protein level in melanoma cell lines. *Acta Biochim. Pol.* 53: 179-188.
6. Chen, H.W., Koehler, C.M. and Teitell, M.A. 2007. Human polynucleotide phosphorylase: location matters. *Trends Cell Biol.* 17: 600-608.
7. Portnoy, V., Palnizky, G., Yehudai-Resheff, S., Glaser, F. and Schuster, G. 2008. Analysis of the human polynucleotide phosphorylase (PNPase) reveals differences in RNA binding and response to phosphate compared to its bacterial and chloroplast counterparts. *RNA* 14: 297-309.
8. Slomovic, S. and Schuster, G. 2008. Stable PNPase RNAi silencing: its effect on the processing and adenylation of human mitochondrial RNA. *RNA* 14: 310-323.

CHROMOSOMAL LOCATION

Genetic locus: PNPT1 (human) mapping to 2p16.1.

PRODUCT

PNPase (h): 293T Lysate represents a lysate of human PNPase transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

PNPase (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive PNPase antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.