

Aph-1 (m2): 293T Lysate: sc-118474

BACKGROUND

Four proteins comprise the γ -secretase complex: Presenilin, nicastrin, Aph-1 and PEN-2. Together, these proteins mediate cell surface signaling pathways for a variety of type I membrane proteins, notably β -Amyloid precursor protein, a protein implicated in the development of Alzheimer's disease, via intramembrane proteolysis. The proteins assemble into a proteolytically active complex in the Golgi/*trans*-Golgi network (TGN) compartments. Assembly leads to autocleavage of Presenilin into two subunits to create the active site of γ -secretase, an important step in understanding the mechanisms involved in the etiology and possible treatment of Alzheimer's disease.

REFERENCES

1. Kimberly, W.T. and Wolfe, M.S. 2003. Identity and function of γ -secretase. *J. Neurosci. Res.* 74: 353-360.
2. Baulac, S., LaVoie, M.J., Kimberly, W.T., Strahle, J., Wolfe, M.S., Selkoe, D.J. and Xia, W. 2003. Functional γ -secretase complex assembly in Golgi/*trans*-Golgi network: interactions among Presenilin, nicastrin, Aph-1, PEN-2, and γ -secretase substrates. *Neurobiol. Dis.* 14: 194-204.
3. Wolfe, M.S. 2003. γ -secretase—intramembrane protease with a complex. *Sci. Aging Knowledge Environ.* 2003: PE7.
4. Fortna, R.R., Crystal, A.S., Morais, V.A., Pijak, D.S., Lee, V.M. and Doms, R.W. 2004. Membrane topology and nicastrin-enhanced endoproteolysis of Aph-1, a component of the γ -secretase complex. *J. Biol. Chem.* 279: 3685-3693.
5. Shiotani, K., Edbauer, D., Prokop, S., Haass, C. and Steiner, H. 2004. Identification of distinct γ -secretase complexes with different Aph-1 variants. *J. Biol. Chem.* 279: 41340-41345.
6. Hansson, E.M., Strömberg, K., Bergstedt, S., Yu, G., Näslund, J., Lundkvist, J. and Lendahl, U. 2005. Aph-1 interacts at the cell surface with proteins in the active γ -secretase complex and membrane-tethered Notch. *J. Neurochem.* 92: 1010-1020.
7. Ma, G., Li, T., Price, D.L. and Wong, P.C. 2005. Aph-1a is the principal mammalian Aph-1 isoform present in γ -secretase complexes during embryonic development. *J. Neurosci.* 25: 192-198.
8. Saito, S. and Araki, W. 2005. Expression profiles of two human Aha-1 genes and their roles in formation of Presenilin complexes. *Biochem. Biophys. Res. Commun.* 327: 18-22.

CHROMOSOMAL LOCATION

Genetic locus: Aph1a (mouse) mapping to 3 F2.1.

PRODUCT

Aph-1 (m2): 293T Lysate represents a lysate of mouse Aph-1 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

Aph-1 (m2): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive Aph-1 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.