



Arg2 (m): 293T Lysate: sc-118522

BACKGROUND

Arginase I (also designated liver-type arginase), which is expressed almost exclusively in the liver, catalyzes the conversion of arginine to ornithine and urea. The human arginase I gene, which maps to chromosome 6q23, encodes a 322 amino acid protein. Arginase I exists as a homotrimeric protein and contains a binuclear manganese cluster. Arginase II catalyzes the same reaction as arginase I, but differs in its tissue specificity and subcellular location. Specifically, arginase II localizes to the mitochondria. Arginase II is expressed in non-hepatic tissues, with the highest levels of expression in the kidneys, but, unlike arginase I, is not expressed in liver. The human arginase II gene, which maps to chromosome 14q24.1, encodes a 354 amino acid protein. In addition, arginase II contains a putative amino-terminal mitochondrial localization sequence.

REFERENCES

1. Diez, A., et al. 1994. Immunological identity of the two different molecular mass constitutive subunits of liver arginase. *Biol. Chem. Hoppe Seyler* 375: 537-541.
2. Gotoh, T., et al. 1996. Molecular cloning of cDNA for nonhepatic mitochondrial arginase (arginase II) and comparison of its induction with nitric oxide synthase in a murine macrophage-like cell line. *FEBS Lett.* 395: 119-122.
3. Gotoh, T., et al. 1997. Chromosomal localization of the human arginase II gene and tissue distribution of its mRNA. *Biochem. Biophys. Res. Commun.* 233: 487-491.
4. Carraway, M.S., et al. 1998. Differential expression of arginase and iNOS in the lung in sepsis. *Exp. Lung Res.* 24: 253-268.
5. Mora, A., et al. 2000. Implications of the S-shaped domain in the quaternary structure of human arginase. *Biochim. Biophys. Acta* 1476: 181-190.
6. Ash, D.E. 2004. Structure and function of arginases. *J. Nutr.* 134: 2760S-2767S.
7. Crombez, E.A. and Cederbaum, S.D. 2004. Hyperargininemia due to liver arginase deficiency. *Mol. Genet. Metab.* 84: 243-251.
8. Ensunsa, J.L., et al. 2004. Reducing arginase activity via dietary manganese deficiency enhances endothelium-dependent vasorelaxation of rat aorta. *Exp. Biol. Med.* 229: 1143-1153.
9. Grodzicki, M., et al. 2004. Arginase activity concentration marking in monitoring of hepatocytes function after orthotopic liver transplantation—preliminary report. *Ann. Transplant.* 9: 54-57.

CHROMOSOMAL LOCATION

Genetic locus: Arg2 (mouse) mapping to 12 C3.

PRODUCT

Arg2 (m): 293T Lysate represents a lysate of mouse Arg2 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

Arg2 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive arginase II antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.