

Atg10 (m): 293T Lysate: sc-118605

BACKGROUND

Autophagy, a process that results in the lysosomal-dependent degradation of cytosolic compartments, is carried out by the autophagosome, which is a double-membrane vesicle whose formation is catalyzed by several autophagy-related gene (Atg) proteins. Atg10 (autophagy-related gene 10), also known as PP12616 or APG10L, is a 220 amino acid protein that localizes to the cytoplasm and plays a role in autophagy, specifically functioning as an E2-like enzyme that provides Atg recognition sites during autophagosome synthesis. Atg10 exists as two isoforms which are produced as a result of alternative splicing events. The gene encoding Atg10 maps to human chromosome 5, which contains 181 million base pairs and comprises nearly 6% of the human genome. Deletion of the p arm of chromosome 5 leads to Cri du chat syndrome, while deletion of the q arm of chromosome 5 altogether is common in therapy-related acute myelogenous leukemias and myelodysplastic syndrome.

REFERENCES

1. Mizushima, N., Yoshimori, T. and Ohsumi, Y. 2002. Mouse APG10 as an APG12-conjugating enzyme: analysis by the conjugation-mediated yeast two-hybrid method. *FEBS Lett.* 532: 450-454.
2. Nemoto, T., Tanida, I., Tanida-Miyake, E., Minematsu-Ikeguchi, N., Yokota, M., Ohsumi, M., Ueno, T. and Kominami, E. 2003. The mouse APG10 homologue, an E2-like enzyme for APG12p conjugation, facilitates MAP-LC3 modification. *J. Biol. Chem.* 278: 39517-39526.
3. Shao, Y., Gao, Z., Feldman, T. and Jiang, X. 2007. Stimulation of Atg12-Atg5 conjugation by ribonucleic acid. *Autophagy* 3: 10-16.
4. Criollo, A., Maiuri, M.C., Tasdemir, E., Vitale, I., Fiebig, A.A., Andrews, D., Molgó, J., Díaz, J., Lavandro, S., Harper, F., Pierron, G., di Stefano, D., Rizzuto, R., Szabadkai, G. and Kroemer, G. 2007. Regulation of autophagy by the inositol trisphosphate receptor. *Cell Death Differ.* 14: 1029-1039.
5. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 610800. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Geng, J. and Klionsky, D.J. 2008. The Atg8 and Atg12 ubiquitin-like conjugation systems in macroautophagy. "Protein modifications: beyond the usual suspects" review series. *EMBO Rep.* 9: 859-864.
7. Williams, R.A., Woods, K.L., Juliano, L., Mottram, J.C. and Coombs, G.H. 2009. Characterization of unusual families of Atg8-like proteins and Atg12 in the protozoan parasite *Leishmania major*. *Autophagy* 5: 159-172.
8. Shin, J.H., Yoshimoto, K., Ohsumi, Y., Jeon, J.S. and An, G. 2009. OsAtg10b, an autophagosome component, is needed for cell survival against oxidative stresses in rice. *Mol. Cells* 27: 67-74.

CHROMOSOMAL LOCATION

Genetic locus: Atg10 (mouse) mapping to 13 C3.

PRODUCT

Atg10 (m): 293T Lysate represents a lysate of mouse P-N transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

Atg10 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive Atg10 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.