Bestrophin-2 (m4): 293T Lysate: sc-118799



The Power to Question

BACKGROUND

Bestrophin-2, also known as BEST2 or VMD2L1 (vitelliform macular dystrophy 2-like protein 1), is a 509 amino acid member of the Bestrophin family of proteins. Members of the Bestrophin family are transmembrane proteins that contain a high percentage of aromatic residues, a conserved RFP (Arg-Phe-Pro) motif and they function as anion channels. Bestrophin-2 forms a calcium-sensitive chloride channel located within the cell membrane. It is also believed that Bestrophin-2 channels may also conduct other physiological anions such as bicarbonate. Bestrophin-2 is mainly expressed in retinal pigment epithelium and colon.

REFERENCES

- Stöhr, H., Marquardt, A., Nanda, I., Schmid, M. and Weber, B.H. 2002. Three novel human VMD2-like genes are members of the evolutionary highly conserved RFP-TM family. Eur. J. Hum. Genet. 10: 281-284.
- Sun, H., Tsunenari, T., Yau, K.W. and Nathans, J. 2002. The vitelliform macular dystrophy protein defines a new family of chloride channels. Proc. Natl. Acad. Sci. USA 99: 4008-4013.
- Tsunenari, T., Sun, H., Williams, J., Cahill, H., Smallwood, P., Yau, K.W. and Nathans, J. 2003. Structure-function analysis of the Bestrophin family of anion channels. J. Biol. Chem. 278: 41114-41125.
- Qu, Z., Fischmeister, R. and Hartzell, C. 2004. Mouse Bestrophin-2 is a bona fide Cl⁻ channel: identification of a residue important in anion binding and conduction. J. Gen. Physiol. 123: 327-340.
- Qu, Z. and Hartzell, C. 2004. Determinants of anion permeation in the second transmembrane domain of the mouse Bestrophin-2 chloride channel. J. Gen. Physiol. 124: 371-382.
- Qu, Z., Chien, L.T., Cui, Y. and Hartzell, H.C. 2006. The anion-selective pore
 of the Bestrophins, a family of chloride channels associated with retinal
 degeneration. J. Neurosci. 26: 5411-5419.
- Pifferi, S., Pascarella, G., Boccaccio, A., Mazzatenta, A., Gustincich, S., Menini, A. and Zucchelli, S. 2006. Bestrophin-2 is a candidate calciumactivated chloride channel involved in olfactory transduction. Proc. Natl. Acad. Sci. USA 103: 12929-12934.

CHROMOSOMAL LOCATION

Genetic locus: Best2 (mouse) mapping to 8 C3.

PRODUCT

Bestrophin-2 (m4): 293T Lysate represents a lysate of mouse Bestrophin-2 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

Bestrophin-2 (m4): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive Bestrophin-2 antibodies. Recommended use: $10\text{-}20~\mu l$ per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

Santa Cruz Biotechnology, Inc. 1.800.457.3801 831.457.3800 fax 831.457.3801 Europe +00800 4573 8000 49 6221 4503 0 www.scbt.com