

PARP-1 (m): 293T Lysate: sc-122382

BACKGROUND

Poly(ADP-ribose) polymerase-1 (PARP-1), also designated PARP, is a nuclear DNA-binding zinc finger protein that influences DNA repair, DNA replication, modulation of chromatin structure and apoptosis. In response to genotoxic stress, PARP-1 catalyzes the transfer of ADP-ribose units from NAD⁺ to a number of acceptor molecules including chromatin. PARP-1 recognizes DNA strand interruptions and can complex with RNA and negatively regulate transcription. Actinomycin D- and etoposide-dependent induction of caspases mediates cleavage of PARP-1 into a p89 fragment that traverses into the cytoplasm. Apoptosis-inducing factor (AIF) translocation from the mitochondria to the nucleus is PARP-1-dependent and is necessary for PARP-1-dependent cell death. PARP-1 deficiencies lead to chromosomal instability due to higher frequencies of chromosome fusions and aneuploidy, suggesting that poly(ADP-ribosyl)ation contributes to the efficient maintenance of genome integrity.

REFERENCES

1. Kaufmann, S.H., et al. 1993. Specific proteolytic cleavage of poly(ADP-ribose) polymerase: an early marker of chemotherapy-induced apoptosis. *Cancer Res.* 53: 2976-3985.
2. Lazebnik, Y.A., et al. 1994. Cleavage of poly(ADP-ribose) polymerase by a proteinase with properties like ICE. *Nature* 371: 346-347.
3. Darmon, A.J., et al. 1995. Activation of the apoptotic protease CPP32 by cytotoxic T cell-derived granzyme B. *Nature* 377: 446-448.
4. Wang, Z.Q., et al. 1997. PARP is important for genomic stability but dispensable in apoptosis. *Genes Dev.* 11: 2347-2358.
5. Jeggo, P.A. 1998. DNA repair: PARP—another guardian angel? *Curr. Biol.* 8: R49-R51.
6. d'Adda di Fagagna, F., et al. 1999. Functions of poly(ADP-ribose) polymerase in controlling telomere length and chromosomal stability. *Nat. Genet.* 23: 76-80.
7. Beneke, R., et al. 2000. DNA excision repair and DNA damage-induced apoptosis are linked to poly(ADP-ribosyl)ation but have different requirements for p53. *Mol. Cell. Biol.* 20: 6695-6703.
8. Vispe, S., et al. 2000. A cellular defense pathway regulating transcription through poly(ADP-ribosyl)ation in response to DNA damage. *Proc. Natl. Acad. Sci. USA* 97: 9886-9891.
9. Herceg, Z. and Wang, Z.Q. 2001. Functions of poly(ADP-ribose) polymerase (PARP) in DNA repair, genomic integrity and cell death. *Mutat. Res.* 477: 97-110.

CHROMOSOMAL LOCATION

Genetic locus: Parp1 (mouse) mapping to 1 H4.

PRODUCT

PARP-1 (m): 293T Lysate represents a lysate of mouse PARP-1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

PARP-1 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive PARP-1 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.