



PHYHIPL (m): 293T Lysate: sc-122551

BACKGROUND

PHYHIPL (phytanoyl-CoA 2-hydroxylase interacting protein-like), also known as phytanoyl-CoA hydroxylase-interacting protein-like, is a 376 amino acid protein that contains one fibronectin type-III domain and belongs to the PHYHIP family. Conserved in chimpanzee, canine, mouse, rat, chicken, zebrafish and *Caenorhabditis elegans*, PHYHIPL exists as three alternatively spliced isoforms. PHYHIPL is a down-regulated target of IRX1, a homeobox tumor suppressor gene linked to gastric carcinoma. PHYHIPL may also play a role in the development of the central system. The gene that encodes PHYHIPL maps to human chromosome 10q21.1.

REFERENCES

1. Gerhard, D.S., Wagner, L., Feingold, E.A., Shenmen, C.M., Grouse, L.H., Schuler, G., Klein, S.L., Old, S., Rasooly, R., Good, P., Guyer, M., Peck, A.M., Derge, J.G., Lipman, D., Collins, F.S., Jang, W., Sherry, S., et al. 2004. The status, quality, and expansion of the NIH full-length cDNA project: the mammalian gene collection (MGC). *Genome Res.* 14: 2121-2127.
2. Deloukas, P., Earthrowl, M.E., Grafham, D.V., Rubenfield, M., French, L., Steward, C.A., Sims, S.K., Jones, M.C., Searle, S., Scott, C., Howe, K., Hunt, S.E., Andrews, T.D., Gilbert, J.G., Swarbreck, D., Ashurst, J.L., et al. 2004. The DNA sequence and comparative analysis of human chromosome 10. *Nature* 429: 375-381.
3. Grupe, A., Li, Y., Rowland, C., Nowotny, P., Hinrichs, A.L., Smemo, S., Kauwe, J.S., Maxwell, T.J., Cherny, S., Doil, L., Tacey, K., van Luchene, R., Myers, A., Wavrant-De Vrièze, F., Kaleem, M., Hollingworth, P., et al. 2006. A scan of chromosome 10 identifies a novel locus showing strong association with late-onset Alzheimer disease. *Am. J. Hum. Genet.* 78: 78-88.
4. Gurok, U., Loebbert, R.W., Meyer, A.H., Mueller, R., Schoemaker, H., Gross, G. and Behl, B. 2007. Laser capture microdissection and microarray analysis of dividing neural progenitor cells from the adult rat hippocampus. *Eur. J. Neurosci.* 26: 1079-1090.
5. Weikard, R., Goldammer, T., Eberlein, A. and Kuehn, C. 2009. Novel transcripts discovered by mining genomic DNA from defined regions of bovine chromosome 6. *BMC Genomics* 10: 186.
6. Guo, X., Liu, W., Pan, Y., Ni, P., Ji, J., Guo, L., Zhang, J., Wu, J., Jiang, J., Chen, X., Cai, Q., Li, J., Zhang, J., Gu, Q., Liu, B., Zhu, Z. and Yu, Y. 2010. Homeobox gene IRX1 is a tumor suppressor gene in gastric carcinoma. *Oncogene* 29: 3908-3920.

CHROMOSOMAL LOCATION

Genetic locus: Phyhipl (mouse) mapping to 10 B5.3.

PRODUCT

PHYHIPL (m): 293T Lysate represents a lysate of mouse PHYHIPL transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

PHYHIPL (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive PHYHIPL antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.