

PRP19 (m): 293T Lysate: sc-122796

BACKGROUND

The spliceosome, the gigantic molecular machine that performs pre-mRNA splicing in eukaryotes, contains over 200 different proteins and 5 RNA molecules (U1, U2, U4, U5 and U6). Pre-mRNA splicing is essential to remove internal non-coding regions of pre-mRNA (introns) and to join the remaining segments (exons) into mRNA before translation. The PRP19-associated complex is required for stable association of U5 and U6 with the spliceosome after U4 is released. Changes within the spliceosome upon binding of the PRP19-associated complex include remodeling of the U6/5' splice site interaction and destabilization of Lsm proteins to allow further interaction of U6 with the intron sequence.

REFERENCES

1. Cheng, S.C., Tarn, W.Y., Tsao, T.Y. and Abelson, J. 1993. PRP19: a novel spliceosomal component. *Mol Cell Biol.* 13: 1876-1882.
2. Tarn, W.Y., Lee, K.R. and Cheng, S.C. 1993. Yeast precursor mRNA processing protein PRP19 associates with the spliceosome concomitant with or just after dissociation of U4 small nuclear RNA. *Proc. Natl. Acad. Sci. USA* 90: 10821-10825.
3. Tarn, W.Y., Hsu, C.H., Huang, K.T., Chen, H.R., Kao, H.Y., Lee, K.R. and Cheng, S.C. 1994. Functional association of essential splicing factors with PRP19 in a protein complex. *EMBO J.* 13: 2421-2431.
4. Nilsen, T.W. 2003. The spliceosome: the most complex macromolecular machine in the cell? *Bioessays* 25: 1147-1149.
5. Turner, I.A., Norman, C.M., Churcher, M.J. and Newman, A.J. 2004. Roles of the U5 snRNP in spliceosome dynamics and catalysis. *Biochem. Soc. Trans.* 32: 928-931.
6. Butcher, S.E. and Brow, D.A. 2005. Towards understanding the catalytic core structure of the spliceosome. *Biochem. Soc. Trans.* 33: 447-449.
7. Valadkhan, S. 2005. snRNAs as the catalysts of pre-mRNA splicing. *Curr. Opin. Chem. Biol.* 9: 603-608.
8. Chan, S.P. and Cheng, S.C. 2005. The Prp19-associated complex is required for specifying and U6 with pre-mRNA during spliceosome activation. *J. Biol. Chem.* 280: 31190-31199.
9. Chen, C.H., Kao, D.I., Chan, S.P., Kao, T.C., Lin, J.Y. and Cheng, S.C. 2006. Functional links between the Prp19-associated complex, U4/U6 biogenesis, and spliceosome recycling. *RNA* 12: 765-774.

CHROMOSOMAL LOCATION

Genetic locus: Prpf19 (mouse) mapping to 19 A.

PRODUCT

PRP19 (m): 293T Lysate represents a lysate of mouse PRP19 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

PRP19 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive PRP19 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.


PRP19 (36-K): sc-101236 is recommended as a positive control antibody for Western Blot analysis of enhanced mouse PRP19 expression in PRP19 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:

1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



PRP19 (36-K): sc-101236. Western blot analysis of PRP19 expression in non-transfected: sc-117752 (A) and mouse PRP19 transfected: sc-122796 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.