



PSPH (m): 293T Lysate: sc-122827

BACKGROUND

Serine is an amino acid required for protein and nucleotide synthesis that may also be involved in cell to cell signaling. PSPH, also known as phosphoserine phosphatase or PSP, is a 225 amino acid Mg^{2+} -dependent enzyme that catalyzes the last and irreversible step in the biosynthesis of serine from carbohydrates, which is the hydrolysis of O-phosphoserine. In the embryonic brain, PSPH is highly expressed in periventricular neural progenitors where it may have a role in neural stem cell proliferation. A lack of PSPH in humans has been shown to cause pre- and postnatal growth retardation as well as moderate psychomotor retardation.

REFERENCES

1. Koch, G.A., Eddy, R.L., Haley, L.L., Byers, M.G., McAvoy, M. and Shows, T.B. 1983. Assignment of the human phosphoserine phosphatase gene (PSP) to the pter leads to q22 region of chromosome 7. *Cytogenet. Cell Genet.* 35: 67-69.
2. Sparkes, R.S., Mohandas, T. and Sparkes, M.C. 1983. The human phosphoserine phosphatase gene (PSP) is mapped to chromosome 7 by somatic cell genetic analysis. *Cytogenet. Cell Genet.* 35: 70-71.
3. Collet, J.F., Gerin, I., Rider, M.H., Veiga-da-Cunha, M. and Van Schaftingen, E. 1997. Human L-3-phosphoserine phosphatase: sequence, expression and evidence for a phosphoenzyme intermediate. *FEBS Lett.* 408: 281-284.
4. Jaeken, J., Detheux, M., Fryns, J.P., Collet, J.F., Alliet, P. and Van Schaftingen, E. 1997. Phosphoserine phosphatase deficiency in a patient with Williams syndrome. *J. Med. Genet.* 34: 594-596.
5. Collet, J.F., Stroobant, V. and Van Schaftingen, E. 1999. Mechanistic studies of phosphoserine phosphatase, an enzyme related to P-type ATPases. *J. Biol. Chem.* 274: 33985-33990.
6. Peeraer, Y., Rabijns, A., Collet, J.F., Van Schaftingen, E. and De Ranter, C. 2004. How calcium inhibits the magnesium-dependent enzyme human phosphoserine phosphatase. *Eur. J. Biochem.* 271: 3421-3427.
7. Tribble, G.D., Mao, S., James, C.E. and Lamont, R.J. 2006. A *Porphyromonas gingivalis* haloacid dehalogenase family phosphatase interacts with human phosphoproteins and is important for invasion. *Proc. Natl. Acad. Sci. USA* 103: 11027-11032.
8. Wang, Y., Maciejewski, B.S., Weissmann, G., Silbert, O., Han, H. and Sanchez-Esteban, J. 2006. DNA microarray reveals novel genes induced by mechanical forces in fetal lung type II epithelial cells. *Pediatr. Res.* 60: 118-124.
9. Nakano, I., Dougherty, J.D., Kim, K., Klement, I., Geschwind, D.H. and Kornblum, H.I. 2007. Phosphoserine phosphatase is expressed in the neuralstem cell niche and regulates neural stem and progenitor cell proliferation. *Stem Cells* 8: 1975-1984.

CHROMOSOMAL LOCATION

Genetic locus: PspH (mouse) mapping to 5 G1.3.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

PSPH (m): 293T Lysate represents a lysate of mouse PSPH transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

PSPH (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive PSPH antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.