

c-Jun (m3): 293T Lysate: sc-125071

BACKGROUND

Genes belonging to the Jun and Fos oncogene families encode nuclear proteins that are found to be associated with a number of transcriptional complexes. The c-Jun protein is a major component of the transcription factor AP-1, originally shown to mediate phorbol ester tumor promoter (TPA)-induced expression of responsive genes through the TPA-response element (TRE). The Jun proteins form homo- and heterodimers which bind the TRE, while Fos proteins are active only as heterodimers with any of the Jun proteins. Fos/Jun heterodimers have a much higher affinity for the TRE than Jun homodimers. Ha-Ras augments c-Jun activity and stimulates phosphorylation of its activation domain. An inhibitor of Fos/Jun function, termed IP-1, associates with Fos and Jun and is inactivated upon phosphorylation induced by the cAMP-dependent protein kinase A (PKA).

REFERENCES

1. Sambucetti, L.C. and Curran, T. 1986. The Fos protein complex is associated with DNA in isolated nuclei and binds to DNA cellulose. *Science* 234: 1417-1419.
2. Bohmann, D., Bos, T.J., Admon, A., Nishimura, T., Vogt, P.K. and Tjian, R. 1987. Human proto-oncogene c-Jun encodes a DNA binding protein with structural and functional properties of transcription factor AP-1. *Science* 238: 1386-1392.
3. Distel, R.J., Ro, H.S., Rosen, B.S., Groves, D.L. and Spiegelman, B.M. 1987. Nucleoprotein complexes that regulate gene expression in adipocyte differentiation: direct participation of c-Fos. *Cell* 49: 835-844.
4. Renz, M., Verrier, B., Kurz, C. and Muller, R. 1987. Chromatin association and DNA-binding properties of the c-Fos proto-oncogene product. *Nucleic Acids Res.* 15: 277-292.
5. Angel, P., Allegretto, E.A., Okino, S.T., Hattori, K., Boyle, W.J., Hunter, T. and Karin, M. 1988. Oncogene Jun encodes a sequence-specific transactivator similar to AP-1. *Nature* 332: 166-171.
6. Franza, B.R., Rauscher, F.J., Josephs, S.F. and Curran, T. 1988. The Fos complex and Fos related antigens recognize sequence elements that contain AP-1 binding sites. *Science* 239: 1150-1153.
7. Auwerx, J. and Sassone-Corsi, P.I. 1991. IP-1: a dominant inhibitor of Fos/Jun whose activity is modulated by phosphorylation. *Cell* 64: 983-993.
8. Binétry, B., Smeal, T. and Karin, M. 1991. Ha-Ras augments c-Jun activity and stimulates phosphorylation of its activation domain. *Nature* 351: 122-127.
9. Waldron, R.T., Whitelegge, J.P., Faull, K.F. and Rozengurt E. 2007. Identification of a novel phosphorylation site in c-Jun directly targeted *in vitro* by protein kinase D. *Biochem. Biophys. Res. Commun.* 356: 361-367.

CHROMOSOMAL LOCATION

Genetic locus: Jun (mouse) mapping to 4 C5.

PRODUCT

c-Jun (m3): 293T Lysate represents a lysate of mouse c-Jun transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

c-Jun (m3): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive c-Jun antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.