

# CyPA (m2): 293T Lysate: sc-125206

## BACKGROUND

The immunosuppressant cyclosporin A (CsA) forms a trimolecular complex with cyclophilin and calcineurins to inhibit calcineurin phosphatase activity. Cyclophilins are conserved, ubiquitous and abundant cytosolic peptidyl-prolyl *cis-trans* isomerases that accelerate the isomerization of XaaPro peptide bonds and the refolding of proteins. Human cyclophilin A (CyPA), an intracellular protein of 165 amino acids, is the target of the CsA and is encoded by a single unique gene conserved from yeast to humans. CyPA is known for its involvement in T cell differentiation and proliferation and is highly expressed in brain. CyPA is incorporated into the virion of the type 1 human immunodeficiency virus (HIV-1) via a direct interaction with the capsid domain of the viral Gag polyprotein and is crucial for efficient viral replication. Cyclophilin B (CyPB) is a member of the cyclophilin family with specific N- and C-terminal extensions. Unlike CyPA, CyPB has a signal sequence leading to its translocation in the endoplasmic reticulum. CyPB is secreted in biological fluids such as blood or milk and binds to a specific receptor present on the human lymphoblastic cell line Jurkat and on human peripheral blood lymphocytes.

## REFERENCES

1. Hasel, K.W., Glass, J.R., Godbout, M. and Sutcliffe, J.G. 1991. An endoplasmic reticulum-specific cyclophilin. *Mol. Cell. Biol.* 11: 3484-3491.
2. Arber, S., Krause, K.H. and Caroni, P. 1992. s-cyclophilin is retained intracellularly via a unique COOH-terminal sequence and colocalizes with the calcium storage protein calreticulin. *J. Cell Biol.* 116: 113-125.
3. Pflugl, G., Kallen, J., Schirmer, T., Jansonius, J.N., Zurini, M.G. and Walkinshaw, M.D. 1993. X-ray structure of a decameric cyclophilin-cyclosporin crystal complex. *Nature* 361: 91-94.
4. Le Hir, M., Su, Q., Weber, L., Woerly, G., Granelli-Piperno, A. and Ryffel, B. 1995. *In situ* detection of cyclosporin A: evidence for nuclear localization of cyclosporine and cyclophilins. *Lab. Invest.* 73: 727-733.
5. Mariller, C., Haendler, B., Allain, F., Denys, A. and Spik, G. 1996. Involvement of the N-terminal part of cyclophilin B in the interaction with specific Jurkat T cell binding sites. *Biochem. J.* 317: 571-576.
6. Vajdos, F.F., Yoo, S., Houseweart, M., Sundquist, W.I. and Hill, C.P. 1997. Crystal structure of cyclophilin A complexed with a binding site peptide from the HIV-1 capsid protein. *Protein Sci.* 6: 2297-2307.
7. Endrich, M.M., Gehrig, P. and Gehring, H. 1999. Maturation-induced conformational changes of HIV-1 capsid protein and identification of two highly affinity sites for cyclophilins in the C-terminal domain. *J. Biol. Chem.* 274: 5326-5332.
8. Wang, P., Cardenas, M.E., Cox, G.E., Perfect, J.R. and Heitman, J. 2001. Two cyclophilin A homologs with shared and distinct functions important for growth and virulence of *Cryptococcus neoformans*. *EMBO Rep.* 2: 511-518.
9. Nahreini, P., Hovland, A.R., Kumar, B., andreatta, C., Edwards-Prasad, J. and Prasad, K.N. 2001. Effects of altered cyclophilin A expression on growth and differentiation of human and mouse neuronal cells. *Cell Mol. Neurobiol.* 21: 65-79.

## CHROMOSOMAL LOCATION

Genetic locus: Ppia (mouse) mapping to 11 A1.

## PRODUCT

CyPA (m2): 293T Lysate represents a lysate of mouse CyPA transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

## APPLICATIONS

CyPA (m2): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive CyPA antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

## STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.