



## DHX33 (m): 293T Lysate: sc-125247

### BACKGROUND

DEAD-box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp, are putative RNA helicases implicated in several cellular processes involving modifications of RNA secondary structure and ribosome/spliceosome assembly. Based on their distribution patterns, some members of this family may be involved in embryogenesis, spermatogenesis and cellular growth and division. DDX33 (DEAD (Asp-Glu-Ala-His) box polypeptide 33), also known as DHX33, is a 707 amino acid nucleolar protein belonging to the DEAD box helicase family. Containing a helicase ATP-binding domain and a helicase C-terminal domain, DDX33 is encoded by a gene located on human chromosome 17. Chromosome 17 comprises over 2.5% of the human genome and encodes over 1,200 genes. Two isoforms of DDX33 exist due to alternative splicing events.

### REFERENCES

1. Schmid, S.R., et al. 1992. D-E-A-D protein family of putative RNA helicases. *Mol. Microbiol.* 6: 283-291.
2. Will, C.L., et al. 2002. Characterization of novel SF3b and 17S U2 snRNP proteins, including a human Prp5p homologue and an SF3b DEAD-box protein. *EMBO J.* 21: 4978-4988.
3. Abdelhaleem, M., et al. 2003. The human DDX and DHX gene families of putative RNA helicases. *Genomics* 81: 618-622.
4. Xu, Y.Z., et al. 2004. Prp5 bridges U1 and U2 snRNPs and enables stable U2 snRNP association with intron RNA. *EMBO J.* 23: 376-385.
5. Cordin, O., et al. 2004. The newly discovered Q motif of DEAD-box RNA helicases regulates RNA-binding and helicase activity. *EMBO J.* 23: 2478-2487.
6. Andersen, J.S., et al. 2005. Nucleolar proteome dynamics. *Nature* 433: 77-83.
7. Linder, P. 2006. Dead-box proteins: a family affair—active and passive players in RNP-remodeling. *Nucleic Acids Res.* 34: 4168-4180.

### CHROMOSOMAL LOCATION

Genetic locus: Dhx33 (mouse) mapping to 11 B4.

### PRODUCT

DHX33 (m): 293T Lysate represents a lysate of mouse DHX33 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

### APPLICATIONS

DHX33 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive DHX33 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

### STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.