

p55 CDC (m): 293T Lysate: sc-125767

BACKGROUND

Cyclins, regulatory subunits which associate with kinases, control many of the important steps in cell cycle progression. The Cdc2 protein kinase (p34Cdc2) exhibits protein kinase activity *in vitro* and exists in a complex with both cyclin B and a protein homologous to p13suc 1. Cdc2 kinase is the active subunit of the M phase promoting factor (MPF) and the M phase-specific Histone H1 kinase. The p34Cdc2/cyclin B complex is required for the G₂ to M transition. An additional cell cycle-dependent protein kinase termed p55 CDC exhibits a high degree of homology with the *S. cerevisiae* proteins Cdc20 and Cdc4. The p55 CDC transcript is readily detectable in a variety of cultured cell lines in growth phase, but disappears when cell growth is chemically arrested. p55 CDC shows kinase activity towards α -casein and myelin basic protein.

REFERENCES

1. Brizuela, L., et al. 1987. p13suc 1 acts in the fission yeast cell division cycle as a component of the p34Cdc2 protein kinase. *EMBO J.* 6: 3507-3514.
2. Arion, D., et al. 1988. Cdc2 is a component of the M phase-specific Histone H1 kinase: evidence for identity with MPF. *Cell* 55: 371-378.
3. Dunphy, W.G., et al. 1988. The *Xenopus* Cdc2 protein is a component of MPF, a cytoplasmic regulator of mitosis. *Cell* 54: 423-431.
4. Morla, A.O., et al. 1989. Reversible tyrosine phosphorylation of Cdc2: dephosphorylation accompanies activation during entry into mitosis. *Cell* 58: 193-203.
5. Pines, J., et al. 1989. Isolation of a human cyclin cDNA: evidence for cyclin mRNA and protein regulation in the cell cycle and for interaction with p34Cdc2. *Cell* 58: 833-846.
6. Jessus, C., et al. 1992. Oscillation of MPF is accompanied by periodic association between Cdc25 and Cdc2-cyclin B. *Cell* 68: 323-332.
7. Weinstein, J., et al. 1994. A novel mammalian protein, p55 CDC, present in dividing cells, is associated with protein kinase activity and has homology to the *Saccharomyces cerevisiae* cell division cycle proteins Cdc20 and Cdc4. *Mol. Cell. Biol.* 14: 3350-3363.
8. Ohtoshi, A., et al. 2000. Human p55 (CDC)/Cdc20 associates with cyclin A and is phosphorylated by the cyclin A-Cdk2 complex. *Biochem. Biophys. Res. Commun.* 268: 530-534.
9. Conway, A.M., et al. 2007. Regulation of neuronal Cdc20 (p55 CDC) expression by the plasticity-related transcription factor zif268. *Synapse* 61: 463-468.

CHROMOSOMAL LOCATION

Genetic locus: Cdc20 (mouse) mapping to 4 D2.1.

PRODUCT

p55 CDC (m): 293T Lysate represents a lysate of mouse p55 CDC transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

p55 CDC (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive p55 CDC antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.