

## PLSCR3 (m): 293T Lysate: sc-125835

### BACKGROUND

The calcium-dependent mitochondrial membrane protein PLSCR3 (phospholipid scramblase 3), also known as PLS3, is a member of the phospholipid scramblase (PLS) family. The PLS family consists of membrane-bound enzymes that participate in the bi-directional movement of phospholipids. PLSCR3 is expressed in a wide variety of tissues but is not found in brain, liver or testis. It is involved in the regulation of mitochondrial respiratory function, morphology and apoptotic responses. More specifically, PLSCR3 is responsible for mediating the transport of cardiolipin from the inner mitochondrial membrane to the outer mitochondrial membrane. Cardiolipin is a major polyglycerophospholipid that plays a role in the regulation of mitochondrial enzymes involved in the generation of ATP. PLSCR3 activity is activated/phosphorylated by PKC  $\delta$  and PLSCR3 functions as a downstream effector in PKC  $\delta$ -induced apoptosis.

### REFERENCES

1. Wiedmer, T., Zhou, Q., Kwok, D.Y. and Sims, P.J. 2000. Identification of three new members of the phospholipid scramblase gene family. *Biochim. Biophys. Acta* 1467: 244-253.
2. Liu, J., Dai, Q., Chen, J., Durrant, D., Freeman, A., Liu, T., Grossman, D. and Lee, R.M. 2003. Phospholipid scramblase 3 controls mitochondrial structure, function, and apoptotic response. *Mol. Cancer Res.* 1: 892-902.
3. Liu, J., Chen, J., Dai, Q. and Lee, R.M. 2003. Phospholipid scramblase 3 is the mitochondrial target of protein kinase C  $\delta$ -induced apoptosis. *Cancer Res.* 63: 1153-1156.
4. Wiedmer, T., Zhao, J., Li, L., Zhou, Q., Hevener, A., Olefsky, J.M., Curtiss, L.K. and Sims, P.J. 2004. Adiposity, dyslipidemia, and insulin resistance in mice with targeted deletion of phospholipid scramblase 3 (PLSCR3). *Proc. Natl. Acad. Sci. USA* 101: 13296-13301.
5. He, Y., Liu, J., Durrant, D., Yang, H.S., Sweatman, T., Lothstein, L. and Lee, R.M. 2005. N-benzyladriamycin-14-valerate (AD198) induces apoptosis through protein kinase C- $\delta$ -induced phosphorylation of phospholipid scramblase 3. *Cancer Res.* 65: 10016-10023.
6. Phillippe, M., Bradley, D.F., Ji, H., Oppenheimer, K.H. and Chien, E.K. 2006. Phospholipid scramblase isoform expression in pregnant rat uterus. *J. Soc. Gynecol. Investig.* 13: 497-501.
7. Van, Q., Liu, J., Lu, B., Feingold, K.R., Shi, Y., Lee, R.M. and Hatch, G.M. 2007. Phospholipid scramblase-3 regulates cardiolipin *de novo* biosynthesis and its resynthesis in growing HeLa cells. *Biochem. J.* 401: 103-109.
8. Mutch, D.M., O'Maille, G., Wikoff, W.R., Wiedmer, T., Sims, P.J. and Siuzdak, G. 2007. Mobilization of pro-inflammatory lipids in obese PLSCR3-deficient mice. *Genome Biol.* 8: R38.
9. He, Y., Liu, J., Grossman, D., Durrant, D., Sweatman, T., Lothstein, L., Epand, R.F., Epand, R.M. and Lee, R.M. 2007. Phosphorylation of mitochondrial phospholipid scramblase 3 by protein kinase C- $\delta$  induces its activation and facilitates mitochondrial targeting of tBid. *J. Cell. Biochem.* 101: 1210-1221.

### CHROMOSOMAL LOCATION

Genetic locus: Plscr3 (mouse) mapping to 11 B3.

### PRODUCT

PLSCR3 (m): 293T Lysate represents a lysate of mouse PLSCR3 transfected 293T cells and is provided as 100  $\mu$ g protein in 200  $\mu$ l SDS-PAGE buffer.

### APPLICATIONS

PLSCR3 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive PLSCR3 antibodies. Recommended use: 10-20  $\mu$ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

### STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.